



# Certificate of Analysis

**Product name**

b-NADH, Reduced, Disodium Salt, Grade II

**Material No.**

10004642103

**Lot No.**

60338834

**Date of Manufacturing**

2022-07-06

**Expiry Date**

2023-07

Parameter	Result	Unit
Appearance	white to yellowish amorphous powder	
Sol.(50 mg/ml water)	clear, colourless to slightly yellowish	
NADH(enzym.)	86,9	%
NADH(abs.340 nm)	88,2	%
NADH(abs.260 nm)	90,7	%
Na(flame photom.)	6,1	%
Water(KF)	3,1	%
NAD(enzym.)	0,54	%
AMP(enzym.)	0,07	%
A 260/A 340	2,30	
A 260/A 240	1,96	
React.rates(LDH)	100	%

**Issue date** 2022-07-06

Page 1 of 2

Roche Diagnostics GmbH

Sandhofer Straße 116; D#68305 Mannheim; Telefon +49#621#759#0; Telefax +49#621#759#2890

Sitz der Gesellschaft: Mannheim - Registergericht: AG Mannheim HRB 3962 - Geschäftsführung: Claus Haberda;

Clemens Schmid - Aufsichtsratsvorsitzender: Dr. Thomas Schinecker



# Certificate of Analysis

**Product name**

b-NADH, Reduced, Disodium Salt, Grade II

**Material No.**

10004642103

**Lot No.**

60338834

**Date of Manufacturing**

2022-07-06

**Expiry Date**

2023-07

**Parameter****Result****Unit**

Ethanol(GC)

3,0

%

This lot meets our specification

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2022-07-06

Page 2 of 2

Roche Diagnostics GmbH

Sandhofer Straße 116; D#68305 Mannheim; Telefon +49#621#759#0; Telefax +49#621#759#2890

Sitz der Gesellschaft: Mannheim - Registergericht: AG Mannheim HRB 3962 - Geschäftsführung: Claus Haberda;

Clemens Schmid - Aufsichtsratsvorsitzender: Dr. Thomas Schinecker