

CERTIFICATE OF ANALYSIS

Product Name: VENTANA anti-ALK (D5F3) Rabbit Monoclonal Primary Antibody
 Ventana REF: 790-4796
 Roche Number: 06687199001
 Batch Number: E02272


Test	Acceptance Criteria	Result
Functional Staining on qualified positive and negative ALK non-small cell lung carcinoma slides	The test batch of anti-ALK (D5F3) Antibody must score 1 or greater on the positive NSCLC tissue and within 0.5 point of the reference batch	PASS
	The negative tissue case(s) must have less than or equal to 0.5 cytoplasmic blushing	PASS
	The background staining observed on the negative reagent control slide must be 0.5 or less	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiry date: 2016-03-28

Manufacture date: 01-Apr-2014

Lucia Romero



08-APR-2014

Quality Department

Date