



a member of the Roche group

## CERTIFICATE OF ANALYSIS

---

Product Name: PAS Staining Kit  
Ventana REF: 860-014  
Roche Number: 05279291001  
Batch Number: F07594

Test	Acceptance Criteria	Result
Functional Staining on at least two of the following: kidney, small intestine, liver, skin, and aorta tissue	Glycogen should stain bright magenta	PASS
	Nuclei should stain purple	PASS
	Negative Control – cells where only the cell nuclei stain purple	PASS

---

This lot met the specifications.

Expiry Date: 2017/02/28

Manufacture Date : 06/NOV/2015

Lucia Romero

*Lucia Romero* *17-Feb-2016*

**Quality Department**

**Date**

**CERTIFICATE OF ANALYSIS**


---

Product: OptiView DAB IHC Detection Kit  
 Ventana **REF**: 760-700  
 Roche Number: 06396500001  
 Batch Number: F08259


Test	Acceptance Criteria	Result
Functional staining on Tonsil and Breast Cancer multiblock and Reference slide	Specificity toward Ki-67 nuclear protein must be assessed by a qualified reader to determine if staining is appropriate (pass). Slides used for qualification of the tissue block should be used as a reference.	PASS
	Background will only be evaluated on the "no antibody" control slides. Background is acceptable if scored $\leq 0.5$ point by a qualified reader.	PASS
Visual Inspection	Reagent must be clear and free of turbidity.	PASS

This lot met the specifications.

Expiry Date: 2017/10/28

Manufacture Date: 04/JAN/2016

---


Jeremy Thomas
12 · JAN · 2016

**Quality Department**
**Date**



*a member of the Roche group*

## **CERTIFICATE OF ANALYSIS**

---

Product Name: BenchMark Special Stains Liquid Coverslip  
Ventana REF: 860-034  
Roche Number: 06523072001  
Batch Number: F08502

Test	Acceptance Criteria	Result
Visual Inspection	Reagent must be clear, colorless, free of turbidity and foreign material	PASS

---

This lot met the specifications.

Expiry date: 2016/12/28

Manufacture date: 11/DEC/2015

Laughlin Hoskinson

*Laughlin Hoskinson* *17-mar-2016*

**Quality Department**

**Date**