



Certificate of Analysis

Product name

ANTI-PAN KERATIN PRIMARY ANTIBODY, 25 ML

Material No.

05266840001

Lot No.

J28647

Date of Manufacturing

2022-11-11

Expiry Date

2024-11-01

Parameter**Result****Unit**

Final Acceptance Testing

Pass

On normal skin, normal liver, breast cancer, prostate cancer, and colon cancer.

The test batch must score 3.5 or greater on normal skin, normal liver, breast cancer, prostate cancer, and colon cancer tissue.

The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result.

The background intensity must be 0.5 or less.

The intensity of the negative reagent control slide must be 0.5 or less.

Bioburden

Pass

Visible bacterial colonies must be 100 or less CFU/mL

Visual Inspection

Pass

Reagent must be clear, and free of turbidity.

Issue date 2022-11-30

Page 1 of 2



Certificate of Analysis

Product name

ANTI-PAN KERATIN PRIMARY ANTIBODY, 25 ML

Material No.

05266840001

Lot No.

J28647

Date of Manufacturing

2022-11-11

Expiry Date

2024-11-01

Parameter**Result****Unit**

Reagent must be free of foreign material and substantially free of intrinsic particulate matter.

This lot meets our specification

Disposition by

Lupe Villarreal

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2022-11-30

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>