

CERTIFICATE OF ANALYSIS


Product Name: CONFIRM anti-bcl-2 (124) Mouse Monoclonal Primary Antibody
 Ventana [REF]: 790-4464
 Roche Number: 05986826001
 Batch Number: E09035

Test	Acceptance Criteria	Result
Functional Staining on qualified tissue slides	Reading against reference slide must score 3.0 or greater on a scale of 0-4+	PASS
	Background reading must score 0.5 or less on a scale of 0-4+	PASS
	Intensity of the negative control slide must be 0.5 or less on a scale of 0-4+	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiry date: 2016-10-28

Manufacture date: 19-NOV-2014


 Patricia Gomez 20/Nov/2014
 Quality Department Date