



Certificate of Analysis

Product name

CONFIRM anti-Desmin (DE-R-11) PAb

Material No.

05267005001

Lot No.

K09229

Date of Manufacturing

2023-05-22

Expiry Date

2025-05-21

Parameter	Result	Unit
Functional Staining	Pass	
Colon, intestine, or other smooth muscle		
CONFIRM anti-Desmin (DE-R-11) Primary Antibody binds to intermediate filaments located in the cytoplasm of muscle cells.		
For specific staining intensity on a 0 to 4 scale, the test batch and reference batch must score 3.5 pt or greater.		
The background intensity must be 0.5 pt or less for the test batch and reference batch.		
The run is invalid if the reference batch does not exhibit acceptable stain intensity.		
The background intensity of the negative reagent control slide must be 0.5 pt or less.		
Visual Inspection	Pass	
Reagent must be clear of foreign material and particulate matter, and free of turbidity		
Bioburden	Pass	

Issue date 2023-06-08

Page 1 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>



Certificate of Analysis

Product name

CONFIRM anti-Desmin (DE-R-11) PAb

Material No.

05267005001

Lot No.

K09229

Date of Manufacturing

2023-05-22

Expiry Date

2025-05-21

Parameter**Result****Unit**

Visible bacterial colonies must be 100 or less CFU/mL

This lot meets our specification

Disposition by

Ken Eblen

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2023-06-08

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>