



Certificate of Analysis

Product name

AMPLIFICATION KIT, CE

Material No.

05266114001

Lot No.

P03866

Date of Manufacturing

2026-04-10

Expiry Date

2028-03-03

Parameter**Result****Unit**

Final Acceptance Testing

Pass

Tissue stained with anti-CD30 and the Amplification Kit must be assigned a staining intensity value greater than or equal to 0.75 to less than or equal to 2 points greater than the nonamplified tissue.

Background should be within less than or equal to 0.5 points of nonamplified tissue.

Tissue stained with anti-CD10 and the Amplification Kit must be assigned a staining intensity value greater than or equal to 0.75 to less than or equal to 2 points greater than the nonamplified tissue.

Background should be within less than or equal to 0.5 points of nonamplified tissue.

Visual Inspection

Pass

Reagent must be clear, and free of turbidity.

Reagent must be free of foreign material and substantially free of intrinsic particulate matter.

This lot meets our specification

Disposition by

Lupe Villarreal

Issue date 2026-04-22

Page 1 of 2



Certificate of Analysis

Product name

AMPLIFICATION KIT, CE

Material No.

05266114001

Lot No.

P03866

Date of Manufacturing

2026-04-10

Expiry Date

2028-03-03

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2026-04-22

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>