

# **9180 Electrolyte Analyzer**

Addendum 1.0 to the Instructions for Use, version 10.0

## Revisions

This addendum concerns a new Reference Electrode (09969772001) and a new Reference Electrode Housing (09969764001) for the 9180 Electrolyte Analyzer.

Please use the *Electrodes* table on page 4 in this addendum instead of the *Electrodes* table in the “Accessories and Consumables” section in *9180 Electrolyte Analyzer Instructions for Use*, version 10.0.

**Table of contents**

Available accessories and consumables ..... 4

Revision 1: New electrode and housing added..... 5

Available accessories and consumables

Attachment Version	Software Version	Revision Date	Changes
1.0	NA	February 2024	First release.

Table 1 Revision history

Roche recommends that you familiarize yourself with the new or revised content provided in this chapter.

The following illustration shows how this content is presented in this document.

	Parameter				
	Cl <sup>-</sup>	Na <sup>+</sup>	K <sup>+</sup>	Ca <sup>2+</sup>	Li <sup>+</sup>
Chloride Electrode 03110451180	X	+	+	+	+
Sodium Electrode 03110419180	+	X	+	+	+
Potassium Electrode 03110338180	+	+	X	+	+
Calcium Electrode 03110354180	+	+	+	X	+
Lithium Electrode 03110664180	+	+	+	+	X
Reference Electrode <sup>(1)</sup> 03112306180	X	X	X	X	X
Reference Electrode Housing 03112284180	X	X	X	X	X
Reference Electrode 9180 09969772001	X	X	X	X	X
Reference Electrode Housing 9180 09969764001	X	X	X	X	X
Dummy Electrode <sup>(2)</sup> 03110192035	+	+	+	+	+

Table 11-1 Electrode order information

(1) Not CE-marked and not RoHS-compliant.

(2) Replaces an unused electrode.

X	Has to be installed
+	Dummy electrode has to be used for proper filling of the measuring chamber

🖨 If you print these pages, Roche recommends printing them single-sided. In this way, you can easily insert the new or revised content in its appropriate location(s) in the existing document.

## Revision 1: New electrode and housing added

*Revised page numbers* This section revises content on page 176 of *9180 Electrolyte Analyzer Instructions for Use*, version 10.0.

