

CERTIFICATE OF ANALYSIS

Product: ISH Protease 1
Ventana [REF]: 780-4147
Roche Number: 05273315001
Batch Number: F01704

Test	Acceptance Criteria	Result
Spectrophotometer testing	Reference and Beginning and End Sample Mean readings must score in the range of 0.50 and 1.30	PASS
Visual Inspection	Product should be clear and colorless, but it is acceptable for the reagent to contain a white flocculent material of a non-biological origin	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiry date: 2017-01-28

Manufacture date: 31-MAR-2015

Laura A. Jimenez  03-APR-2015
Quality Department Date