



Certificate of Analysis

Product name

Anti-p63 (4A4) Mouse Monoclonal Primary Antibody

Material No.

05867061001

Lot No.

G09114

Date of Manufacturing

28/Nov/2016

Expiry Date

28/Nov/2018

Parameter**Result****Unit**

Visual Inspection

Pass

Reagent must be free of foreign material and substantially free of particulate matter

Reagent must be clear and free of turbidity

Functional Staining

Pass

Prostate carcinoma and tonsil

The test batch must score 3 or greater on prostate carcinoma and tonsil tissue

The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result

The background intensity must be within 0.5 of the negative control slide

The intensity of the negative reagent control slide must be 0.5 or less

Bioburden

Pass

Visible bacterial colonies must be 100 or less CFU/mL

This lot meets our specification

Disposition by

Laughlin Hoskinson

Issue date

Dec 2016