



a member of the Roche group

CERTIFICATE OF ANALYSIS

Product Name: JONES LIGHT GREEN STAINING KIT

Ventana **REF** : 860-020

Roche Number: 05279356001

Batch Number: F01379

Test	Acceptance Criteria	Result
Slide Staining Test	The basement membrane in kidney stains black while all other tissue structures are light green. Fine basement membrane should be present in the glomeruli of the kidney sections. Negative control- Light green background staining acts as a negative contro	PASS
Final Container Test	Staining comparable to the reference slide.	PASS

This lot met the specifications.

Expiry date: 2017-01-28

Manufacture date : 23-Apr-2015

Jeremy Thomas

24/APR/2015

Quality Department

Date