

CERTIFICATE OF ANALYSIS

Product name: anti-CD30 (Ber-H2) Mouse Monoclonal Primary Antibody
Ventana REF : 790-4858
Roche Number: 07007841001
Batch Number: F06578

Test	Acceptance Criteria	Result
Functional Staining on positive Mature T-Cell Lymphoma, Diffuse Large B-cell Lymphoma, and Hodgkin's lymphoma	-- The test batch of anti-CD30 (Ber-H2) must score greater than or equal to 2.5 on the positive tissue and less than or equal to 0.5 on the negative tissue (scale of 0-3+).	Pass
	--. The test antibody must stain within a 0.5 staining intensity when compared to the reference slide.	Pass
	-- The background staining observed on the positive and negative slides must be 0.5 or less.	Pass
Visual Inspection	Reagent must be free of foreign material and substantially free of intrinsic matter.	Pass
	Reagent must be clear and free of turbidity.	Pass
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	Pass

This lot met the specifications.
Expiry Date: 2016/09/28
Date Manufactured: 11-Sep-2015
Guadalupe Villarreal
Quality Department
08 OCT 2015
Date