

## Customer Address:

Customer Number : 2049990013  
 Customer Mat # : 06523897001  
 Roche Diagnostics Corporation  
 Accounts Payable  
 Po Box 50457  
 INDIANAPOLIS IN 46250-0417  
 UNITED STATES

## Certificate of Analysis

**Product Name :** anti-p57/Kip2 (Kp10) Mouse Monoclonal Primary Antibody  
**Product Number :** 760-4617  
**Batch Number :** V0003672  
 Source Batch : 0000207802  
 Storage Temperature : 2°C - 8°C  
 Expiry Date : 2025-10-31  
 Quality Release Date : 2022-12-22  
 Date of Manufacture : 2022-11-30

Test	Specification	Result
Appearance (Turbidity)	Free of Visible Foreign Material	Free of Visible Foreign Material
Appearance (Color)	Green Transparent	Green Transparent
Appearance (Form)	Liquid	Liquid
IHC Staining	Pass	Pass
Intensity difference between sample and a reference must be less than or equal to 0.5.		
IHC Staining	Pass	Pass
Background difference between sample and a reference must be less than or equal to 0.5.		
IHC - Tissue	Tissue Type	Placenta
ELISA Titer	1.0 - 10.0 µg/ml	1.0 µg/ml
Bioburden	Pass	Pass
Visible bacterial colonies must be less than 100 CFU/mL.		
Source	Supernatant	Supernatant

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.



**Product** 760-4617  
**Batch** V0003672



JUNIOR NUNEZ, Quality, QC Lab Technician II  
Quality Assurance  
Rocklin, California  
US

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.

Version Number: 01 Doc: 1096197

Page 2 of 2

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact [mlsbranding@sial.com](mailto:mlsbranding@sial.com)

