



## Certificate of Analysis

**Product name**

KIT DNA ISOLATION 24T IVD

**Material No.**

05985536190

**Lot No.**

E25665

**Date of Manufacturing**

2019-01-26

**Expiry Date**

2020-09

Parameter	Result	Unit
Proc K562 Reps-NMD	9	#lnRng
Proc K562-Ch0 Ct	9	#lnRng
Proc K562-Ch1 Ct	9	#lnRng
Proc K562-Ch1 Avg	25.0	MEAN
DNASPKit-IVD	PASS	

The visual inspection should be as follows:

1 x DNA TLB (DNA Tissue Lysis Buffer, 10 mL)  
Natural bottle and screw cap containing a clear,  
colorless liquid

1 x DNA PBB (DNA Paraffin Binding Buffer, 10 mL)  
Natural bottle and screw cap containing a clear,  
colorless liquid

1 x WB I (FFPET Wash Buffer I, 25 mL)  
Natural bottle and screw cap containing a clear,  
colorless liquid

1 x WB II (FFPET Wash Buffer II, 12.5 mL)

**Issue date** 2019-01-31

Page 1 of 3

Roche Molecular Systems, Inc.  
Quality Management  
1080 U.S. Highway 202 South  
Branchburg, NJ 08876-3771



## Certificate of Analysis

**Product name**

KIT DNA ISOLATION 24T IVD

**Material No.**

05985536190

**Lot No.**

E25665

**Date of Manufacturing**

2019-01-26

**Expiry Date**

2020-09

**Parameter****Result****Unit**

Natural bottle and screw cap containing a clear, colorless liquid

1 x DNA EB (DNA Elution Buffer, 6 mL)

Natural bottle and screw cap containing a clear, colorless liquid

1 x PK (Proteinase K, 100 mg)

Brown glass bottle containing a white, clumpy powder

1 x FT (Filter Tubes with Caps, 25 pcs)

Bag containing 25 filter tubes with caps

3 x CT (Collection Tubes, 25 pcs)

Bag containing 25 polypropylene tubes

Components

1

This lot meets our specification

**Disposition by**

Craig Cawley

**Issue date** 2019-01-31

Page 2 of 3

Roche Molecular Systems, Inc.  
Quality Management  
1080 U.S. Highway 202 South  
Branchburg, NJ 08876-3771



## Certificate of Analysis

**Product name**

KIT DNA ISOLATION 24T IVD

**Material No.**

05985536190

**Lot No.**

E25665

**Date of Manufacturing**

2019-01-26

**Expiry Date**

2020-09

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2019-01-31

Page 3 of 3

Roche Molecular Systems, Inc.  
Quality Management  
1080 U.S. Highway 202 South  
Branchburg, NJ 08876-3771