

## Customer Address:

Customer Number : 2049990013  
 Customer Mat # : 05269784001  
 Roche Diagnostics Corporation  
 Accounts Payable  
 Po Box 50457  
 INDIANAPOLIS IN 46250-0417  
 UNITED STATES

## Certificate of Analysis


**Product Name :** Glial Fibrillary Acidic Protein (EP672Y) Rabbit Monoclonal Antibody  
**Product Number :** 760-4345  
**Batch Number :** V0004200  
 Source Batch : 0000269615  
 Storage Temperature : 2°C - 8°C  
 Expiry Date : 2026-08-31  
 Quality Release Date : 2023-11-17  
 Date of Manufacture : 2023-09-22

Test	Specification	Result
Appearance (Turbidity)	Free of Visible Foreign Material	Free of Visible Foreign Material
Appearance (Color)	Green Transparent	Green Transparent
Appearance (Form)	Liquid	Liquid
IHC Staining Intensity difference between sample and a reference must be less than or equal to 0.5.	Pass	Pass
IHC Staining Background difference between sample and a reference must be less than or equal to 0.5.	Pass	Pass
IHC - Tissue	Tissue Type	Brain
ELISA Titer	1 - 10 µg/ml	1 µg/ml
Bioburden Visible bacterial colonies must be less than 100 CFU/mL.	Pass	Pass
Source	Supernatant	Supernatant

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.



**Product** 760-4345  
**Batch** V0004200



Bo Vang, Quality, QC Lab Technician I  
Quality Assurance  
Rocklin, California  
US

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.

