

Certificate of Analysis

Product Name: GH (Growth Hormone) Rabbit Polyclonal Antibody
Catalog Number(s): 760-2804
Roche Part Number: 05268257001
Lot Number: 1613704B
Date of Manufacture: 2016-05-16
Expiration Date: 2019-04-30

Test	Acceptance Criteria	Results
Functional testing on the following tissue(s): Pituitary	Intensity difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
	Background difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Appearance	Transparent liquid and free of visible foreign material.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Bioburden	Visible bacterial colonies must be less than 100CFU/mL.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
ELISA	2.84-5.28 µg/mL	4.08 µg/mL
Source	N/A	Antisera

Additional information: N/A



CELL MARQUE
6600 Sierra College Blvd.
Rocklin, California 95677
(916) 746-8900 Fax: (916) 746-8989

Authorized Signature:  Date: 08/10/17
Quality Assurance Representative

Issued: 6/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]