

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**SECTION 1. IDENTIFICATION**

Product name : CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary Antibody

Product code : 05278384001

Manufacturer or supplier's details

Company name of supplier : Roche Diagnostics GmbH

Address : 9115 Hague Road
Indianapolis, IN 46250
USA

Telephone : 1-800-428-5074

Emergency telephone

In case of emergencies : CHEMTREC 1-800-424-9300 (U.S. or Canada)
1-703-527-3887 (International)**Recommended use of the chemical and restrictions on use**Recommended use : Laboratory chemicals
Refer to product literature for further details.**SECTION 2. HAZARDS IDENTIFICATION****GHS classification in accordance with the OSHA Hazard Communication Standard (29 CFR 1910.1200)**

Skin sensitization : Category 1

GHS label elements

Hazard pictograms :



Signal Word : Warning

Hazard Statements : H317 May cause an allergic skin reaction.

Precautionary Statements : **Prevention:**
P261 Avoid breathing mist or vapors.
P272 Contaminated work clothing must not be allowed out of the workplace.
P280 Wear protective gloves.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**Response:**

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.
P333 + P313 If skin irritation or rash occurs: Get medical advice/
attention.
P363 Wash contaminated clothing before reuse.

Disposal:

P501 Dispose of contents/ container to an approved waste disposal plant.

Other hazards

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Components

Chemical name	CAS-No.	Concentration (% w/w)
Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H - isothiazol-3-one (3:1)	55965-84-9	≥ 0.0015 - < 0.06
L-Ascorbic acid	50-81-7	< 0.1
L-Tryptophan	73-22-3	< 0.1
Vitamin B12	68-19-9	< 0.1
L-Ornithine, hydrochloride (1:1)	3184-13-2	< 0.1
Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1)	67-48-1	< 0.1

Actual concentration is withheld as a trade secret

SECTION 4. FIRST AID MEASURES

General advice : Move out of dangerous area.
Show this material safety data sheet to the doctor in attendance.
Do not leave the victim unattended.

If inhaled : If unconscious, place in recovery position and seek medical advice.
If symptoms persist, call a physician.

In case of skin contact : If on skin, rinse well with water.

In case of eye contact : Flush eyes with water as a precaution.
Remove contact lenses.
Protect unharmed eye.
Keep eye wide open while rinsing.
If eye irritation persists, consult a specialist.

If swallowed : Keep respiratory tract clear.
Do not give milk or alcoholic beverages.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Never give anything by mouth to an unconscious person.
If symptoms persist, call a physician.

Most important symptoms and effects, both acute and delayed : May cause an allergic skin reaction.

Notes to physician : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media : High volume water jet

Specific hazards during fire fighting : Do not allow run-off from fire fighting to enter drains or water courses.

Hazardous combustion products : Carbon oxides

Further information : Collect contaminated fire extinguishing water separately. This must not be discharged into drains.
Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations.

Special protective equipment for fire-fighters : Wear self-contained breathing apparatus for firefighting if necessary.

SECTION 6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures : Use personal protective equipment.
Refer to protective measures listed in sections 7 and 8.

Environmental precautions : Prevent product from entering drains.
Prevent further leakage or spillage if safe to do so.
If the product contaminates rivers and lakes or drains inform respective authorities.

Methods and materials for containment and cleaning up : Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder, sawdust).



CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary Anti-body

Version
2.1

Revision Date:
02-25-2025

Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

- Advice on protection against fire and explosion : Normal measures for preventive fire protection.
- Advice on safe handling : Do not breathe vapors/dust.
Avoid exposure - obtain special instructions before use.
Avoid contact with skin and eyes.
For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.
Dispose of rinse water in accordance with local and national regulations.
Persons susceptible to skin sensitization problems or asthma, allergies, chronic or recurrent respiratory disease should not be employed in any process in which this mixture is being used.
- Conditions for safe storage : Keep container tightly closed in a dry and well-ventilated place.
Containers which are opened must be carefully resealed and kept upright to prevent leakage.
Electrical installations / working materials must comply with the technological safety standards.
- Further information on storage conditions : See label, package insert or internal guidelines
- Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
L-Ascorbic acid	50-81-7	IOEL	10 mg/m3	Roche Industrial Hygiene Committee (RIHC)
L-Tryptophan	73-22-3	IOEL	10 mg/m3	Roche Industrial Hygiene Committee (RIHC)
Vitamin B12	68-19-9	IOEL	0.6 mg/m3	Roche In-



CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary Anti-body

Version
2.1

Revision Date:
02-25-2025

Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

				dustrial Hy- giene Com- mittee (RIHC)
L-Ornithine, hydrochloride (1:1)	3184-13-2	IOEL	10 mg/m3	Roche In- dustrial Hy- giene Com- mittee (RIHC)
		IOEL	10 mg/m3	Roche Indus- trial Hygiene Committee (RIHC)
Ethanaminium, 2-hydroxy- N,N,N-trimethyl-, chloride (1:1)	67-48-1	IOEL	10 mg/m3	Roche In- dustrial Hy- giene Com- mittee (RIHC)

Engineering measures : No data available

Personal protective equipment

Respiratory protection : General and local exhaust ventilation is recommended to maintain vapor exposures below recommended limits. Where concentrations are above recommended limits or are unknown, appropriate respiratory protection should be worn. Follow OSHA respirator regulations (29 CFR 1910.134) and use NIOSH/MSHA approved respirators. Protection provided by air purifying respirators against exposure to any hazardous chemical is limited. Use a positive pressure air supplied respirator if there is any potential for uncontrolled release, exposure levels are unknown, or any other circumstance where air purifying respirators may not provide adequate protection.

Hand protection

In case of contact through splashing:
Material : Nitrile rubber
Break through time : > 30 min
Glove thickness : > 0.11 mm

In case of full contact:
Material : butyl-rubber
Break through time : > 480 min
Glove thickness : > 0.4 mm

Remarks : Wear appropriate protective gloves to prevent skin contact.
Replace torn or punctured gloves promptly.

Eye protection : Eye wash bottle with pure water
Tightly fitting safety goggles

Skin and body protection : Impervious clothing

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Choose body protection according to the amount and concentration of the dangerous substance at the work place.

Hygiene measures : Wash hands before breaks and at the end of workday.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	: liquid
Color	: colorless
Odor	: odorless
Odor Threshold	: No data available
pH	: 7.7 Concentration: 100 %
Melting point/ range	: No data available
Boiling point/boiling range	: 212 °F / 100 °C
Flash point	: does not flash
Evaporation rate	: No data available
Flammability (liquids)	: Does not sustain combustion.
Self-ignition	: Not applicable
Upper explosion limit / Upper flammability limit	: No data available
Lower explosion limit / Lower flammability limit	: No data available
Vapor pressure	: No data available
Relative vapor density	: No data available
Relative density	: No data available
Density	: 1.005 g/cm ³ (68 °F / 20 °C, 1,013 hPa)
Solubility(ies) Water solubility	: completely miscible

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Solubility in other solvents	:	No data available
Partition coefficient: n-octanol/water	:	No data available
Autoignition temperature	:	No data available
Decomposition temperature	:	No data available
Viscosity		
Viscosity, dynamic	:	No data available
Viscosity, kinematic	:	No data available
Explosive properties	:	Not explosive
Oxidizing properties	:	The substance or mixture is not classified as oxidizing.
Particle characteristics		
Particle Size Distribution	:	Not applicable

SECTION 10. STABILITY AND REACTIVITY

Reactivity	:	No dangerous reaction known under conditions of normal use.
Chemical stability	:	Stable under normal conditions.
Possibility of hazardous reactions	:	No decomposition if stored and applied as directed.
Conditions to avoid	:	No data available
Incompatible materials	:	Not applicable
Hazardous decomposition products	:	No hazardous decomposition products are known.

SECTION 11. TOXICOLOGICAL INFORMATION**Acute toxicity**

Not classified due to lack of data.

Product:

Acute inhalation toxicity	:	Acute toxicity estimate: > 200 mg/l Exposure time: 4 h Test atmosphere: vapor Method: Calculation method
---------------------------	---	-------------------------------------------------------------------------------------------------------------------

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**Components:****Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1):**

Acute oral toxicity : LD50 Oral (Rat, male): 64 mg/kg

Acute inhalation toxicity : LC50 (Rat, male and female): 0.33 mg/l
Exposure time: 4 h
Test atmosphere: dust/mist
Method: OECD Test Guideline 403

Acute dermal toxicity : LD50 Dermal (Rabbit, male): 87.12 mg/kg

L-Ascorbic acid:

Acute oral toxicity : LD50 Oral (Rat): 11,900 mg/kg
LD50 Oral (Mouse): 8,000 mg/kg

L-Tryptophan:

Acute oral toxicity : LD50 Oral (Rat): > 16,000 mg/kg

Vitamin B12:

Acute oral toxicity : LD50 Oral (Mouse): > 5,000 mg/kg

L-Ornithine, hydrochloride (1:1):

Acute oral toxicity : LD50 Oral (Rat): 10,000 mg/kg

Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):

Acute oral toxicity : LD50 Oral (Rat): 3,400 mg/kg
LD50 Oral (Mouse): 3,900 mg/kg
Acute toxicity estimate: > 5,000 mg/kg
Method: Expert judgment

Acute inhalation toxicity : Acute toxicity estimate: > 30 mg/l
Test atmosphere: dust/mist

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Method: Expert judgment

Acute dermal toxicity : Acute toxicity estimate: > 5,000 mg/kg
Method: Expert judgment**Skin corrosion/irritation**

Not classified due to lack of data.

Product:

Remarks : May cause skin irritation and/or dermatitis.

Components:**Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1):**

Result : Corrosive after 1 to 4 hours of exposure

L-Ascorbic acid:

Remarks : This information is not available.

L-Tryptophan:Species : Rabbit
Method : OECD Test Guideline 404
Result : No skin irritation**Vitamin B12:**

Remarks : May cause skin irritation in susceptible persons.

Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):

Remarks : This information is not available.

Serious eye damage/eye irritation

Not classified due to lack of data.

Product:

Remarks : Vapors may cause irritation to the eyes, respiratory system and the skin.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**Components:****Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1):**

Species : Rabbit
Result : Risk of serious damage to eyes.

L-Tryptophan:

Species : Rabbit
Result : No eye irritation
Method : OECD Test Guideline 405

Vitamin B12:

Remarks : May irritate eyes.

L-Ornithine, hydrochloride (1:1):

Result : Irritating to eyes.

Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):

Remarks : This information is not available.

Respiratory or skin sensitization**Skin sensitization**

May cause an allergic skin reaction.

Respiratory sensitization

Not classified due to lack of data.

Product:

Remarks : Causes sensitization.

Components:**Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1):**

Test Type : Local lymph node assay (LLNA)
Species : Mouse
Assessment : The product is a skin sensitizer, sub-category 1A.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**L-Tryptophan:**

Test Type : Local lymph node assay (LLNA)
Species : Mouse
Method : OECD Test Guideline 429
Result : Does not cause skin sensitization.

Vitamin B12:

Remarks : May cause sensitization of susceptible persons.

Germ cell mutagenicity

Not classified due to lack of data.

Components:**L-Tryptophan:**

Genotoxicity in vitro : Test Type: reverse mutation assay
Test system: Salmonella typhimurium
Metabolic activation: with and without metabolic activation
Method: OECD Test Guideline 471
Result: negative

Vitamin B12:

Genotoxicity in vitro : Test Type: reverse mutation assay
Metabolic activation: with and without metabolic activation
Method: OECD Test Guideline 471
Result: negative

Carcinogenicity

Not classified due to lack of data.

Components:**L-Ascorbic acid:**

Remarks : No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

L-Tryptophan:

Remarks : No ingredient of this product present at levels greater than or

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

L-Ornithine, hydrochloride (1:1):

Remarks : No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

Reproductive toxicity

Not classified due to lack of data.

STOT-single exposure

Not classified due to lack of data.

Components:**L-Ascorbic acid:**

Assessment : The substance or mixture is not classified as specific target organ toxicant, single exposure.

L-Tryptophan:

Assessment : The substance or mixture is not classified as specific target organ toxicant, single exposure.

Vitamin B12:

Assessment : The substance or mixture is not classified as specific target organ toxicant, single exposure.

STOT-repeated exposure

Not classified due to lack of data.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**Components:****L-Ascorbic acid:**

Assessment : The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

L-Tryptophan:

Assessment : The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Vitamin B12:

Assessment : The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):

Assessment : The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Repeated dose toxicity**Components:****L-Tryptophan:**

Species : Rat, female
NOAEL : ca. 3764 mg/kg
Application Route : Oral
Exposure time : 90 d
Dose : 0, 401, 1214 or 3946
Method : OECD Test Guideline 408
GLP : yes

Aspiration toxicity

Not classified due to lack of data.

Components:**L-Ascorbic acid:**

No data available

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024**L-Tryptophan:**

No data available

Vitamin B12:

No data available

Further information**Product:**

Remarks : No data available

Components:**Vitamin B12:**

Remarks : May be harmful by inhalation, ingestion, skin adsorption.

SECTION 12. ECOLOGICAL INFORMATION**Ecotoxicity****Components:****Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1):**Toxicity to fish : LC50 (Oncorhynchus mykiss (rainbow trout)): 0.19 mg/l
Exposure time: 96 h
Test Type: flow-through testToxicity to daphnia and other : EC50 (Daphnia magna (Water flea)): 0.16 mg/l
aquatic invertebrates Exposure time: 48 h
Test Type: flow-through testToxicity to algae/aquatic : ErC50 (Skeletonema costatum (marine diatom)): 0.0063 mg/l
plants Exposure time: 72 h
Test Type: static test
Method: OECD Test Guideline 201M-Factor (Acute aquatic tox- : 100
icity)Toxicity to fish (Chronic tox- : NOEC (Cyprinodon sp. (minnow)): 0.02 mg/l
icity) Exposure time: 38 d
Test Type: semi-static test
Method: OECD Test Guideline 210Toxicity to daphnia and other : NOEC (Daphnia): 0.1 mg/l
aquatic invertebrates (Chron- Exposure time: 21 d
ic toxicity) GLP: yes

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

M-Factor (Chronic aquatic toxicity) : 100

Ecotoxicology Assessment

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

L-Ascorbic acid:Toxicity to fish : LC50 (Oncorhynchus mykiss (rainbow trout)): 1,020 mg/l
Exposure time: 96 h
Method: OECD Test Guideline 203**Ecotoxicology Assessment**

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

L-Tryptophan:**Ecotoxicology Assessment**

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Vitamin B12:**Ecotoxicology Assessment**

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):Toxicity to fish : LC50 (Fish): > 100 mg/l
Exposure time: 96 h**Ecotoxicology Assessment**

Acute aquatic toxicity : This product has no known ecotoxicological effects.

Chronic aquatic toxicity : This product has no known ecotoxicological effects.

Toxicity Data on Soil : Not expected to adsorb on soil.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Other organisms relevant to the environment : No data available

Persistence and degradability**Components:****Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1):**

Biodegradability : Result: Inherently biodegradable.
Remarks: The 10 day time window criterion is not fulfilled.

Result: Not readily biodegradable.
Method: OECD Test Guideline 301B
GLP: yes

L-Ascorbic acid:

Biodegradability : Biodegradation: 97 %
Exposure time: 5 d
Method: OECD Test Guideline 302B

Biodegradation: 100 %
Exposure time: 150 d
Method: OECD Test Guideline 302B

Vitamin B12:

Biodegradability : Result: Not readily biodegradable.

Bioaccumulative potential**Components:****Mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1):**

Bioaccumulation : Species: Lepomis macrochirus (Bluegill sunfish)
Bioconcentration factor (BCF): 54
Temperature: 68 °F / 20 °C
Concentration: ca. 0.01 mg/l

Partition coefficient: n-octanol/water : log Pow: 0.75

L-Ascorbic acid:

Partition coefficient: n-octanol/water : log Pow: -2.15 (73 °F / 23 °C)

L-Tryptophan:

Partition coefficient: n-octanol/water : log Pow: -1.05

Vitamin B12:

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Bioaccumulation : Species: Fish
Bioconcentration factor (BCF): 3
Exposure time: 3 d

Partition coefficient: n-octanol/water : log Pow: ca. 3.57

L-Ornithine, hydrochloride (1:1):

Partition coefficient: n-octanol/water : Remarks: No data available

Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):

Partition coefficient: n-octanol/water : Remarks: No data available

Mobility in soil

No data available

Other adverse effects**Product:**

Ozone-Depletion Potential : Regulation: 40 CFR Protection of Environment; Part 82 Protection of Stratospheric Ozone - CAA Section 602 Class I Substances
Remarks: This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

Additional ecological information : An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.
Harmful to aquatic life with long lasting effects.

Components:**Ethanaminium, 2-hydroxy-N,N,N-trimethyl-, chloride (1:1):**

Adsorbed organic bound halogens (AOX) : Remarks: Not applicable

Additional ecological information : No data available

SECTION 13. DISPOSAL CONSIDERATIONS**Disposal methods**

Waste from residues : The product should not be allowed to enter drains, water courses or the soil.
Do not contaminate ponds, waterways or ditches with chemical or used container.
Send to a licensed waste management company.

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary Anti-bodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Contaminated packaging : Empty remaining contents.
Dispose of as unused product.
Do not re-use empty containers.

SECTION 14. TRANSPORT INFORMATION**International Regulations****UNRTDG**

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

Domestic regulation**49 CFR**

Not regulated as a dangerous good

Special precautions for user

Remarks : Not dangerous goods in the meaning of ADR/RID, ADN,
IMDG-Code, ICAO/IATA-DGR

SECTION 15. REGULATORY INFORMATION**CERCLA Reportable Quantity**

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 304 Extremely Hazardous Substances Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

This material does not contain any components with a section 302 EHS TPQ.

SARA 311/312 Hazards : Respiratory or skin sensitization

SARA 313 : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Clean Air Act

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).
This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 111 SOCM Intermediate or Final VOC's (40 CFR 60.489).

Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

This product does not contain any priority pollutants related to the U.S. Clean Water Act

US State Regulations**Massachusetts Right To Know**

Sodium azide (Na(N₃)) 26628-22-8

Pennsylvania Right To Know

Water 7732-18-5

Maine Chemicals of High Concern

Product does not contain any listed chemicals

Vermont Chemicals of High Concern

Vitamin B12 68-19-9

Washington Chemicals of High Concern

Vitamin B12 68-19-9

The ingredients of this product are reported in the following inventories:

AIIC : Not in compliance with the inventory

DSL : This product contains the following components that are not on the Canadian DSL nor NDSL.

Animal serum /-plasma - different species

Glycols

Serum, Calf

Modified alkyl carboxylate

L-Histidine monohydrochloride monohydrate

L-Tyrosine disodiumsalt dihydrate

2'-Deoxyadenosine monohydrate

N-methyltaurine

CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary AntibodyVersion
2.1Revision Date:
02-25-2025Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

Uridine-5'-diphospho-glucuronate, Disodium salt (UDP-glucuronate)

Tetramethylrhodamin-6-2'-desoxy-uridin-5-triphosphate

Coenzyme A

D-Glucuronic acid

Riboflavin 5'-(trihydrogen diphosphate), P'.->.5'-ester with adenosine

NZIoC	: Not in compliance with the inventory
ENCS	: Not in compliance with the inventory
ISHL	: Not in compliance with the inventory
KECI	: Not in compliance with the inventory
PICCS	: Not in compliance with the inventory
IECSC	: Not in compliance with the inventory
TCSI	: Not in compliance with the inventory
TSCA	: Product contains substance(s) not listed on TSCA inventory.
TECI	: Not in compliance with the inventory

TSCA list

No substances are subject to a Significant New Use Rule.

No substances are subject to TSCA 12(b) export notification requirements.

SECTION 16. OTHER INFORMATION**Further information**

SAFETY DATA SHEET




CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary Antibody

Version
2.1

Revision Date:
02-25-2025

Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

NFPA 704:



HMIS® IV:

HEALTH	/	2
FLAMMABILITY		0
PHYSICAL HAZARD		0

HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. The "*" represents a chronic hazard, while the "/" represents the absence of a chronic hazard.

Full text of other abbreviations

AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECL - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act

SAFETY DATA SHEET



CONFIRM anti-Ki-67 (30-9) Rabbit Monoclonal Primary Antibody

Version
2.1

Revision Date:
02-25-2025

Date of last issue: 06-28-2024
Date of first issue: 03-19-2024

(United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 02-25-2025

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

US / Z8 / 2404