



Document Number
Form 7.5-58-2

Title
Certificate of Analysis - Predilute

Authoring Site	Document Type – Subtype	Document Security
Rocklin, CA	Procedural – Form	N/A
Referenced Documents		Special Notes
1. WI 7.5-58 Releasing Finished Goods Inventory in SAP		N/A

Certificate of Analysis

Product Name: GATA3 (L50-823) Mouse Monoclonal Primary Antibody
Product Number: 760-4897
Batch Number: V0002352
 Storage Temperature: 2°C - 8°C
 Expiration Date: 2023-07-31
 Quality Release Date: 2020-10-08
 Date of Manufacture: 2020-08-25

Test	Specification	Result
Appearance (Turbidity)	Free of Visible Foreign Material	Free of Visible Foreign Material
Appearance (Color)	Green Transparent	Green Transparent
Appearance (Form)	Liquid	Liquid
IHC Staining	Pass	Pass
Intensity difference between sample and a reference must be less than or equal to 0.5		
IHC Staining	Pass	Pass
Background difference between sample and a reference must be less than or equal to 0.5		
IHC – Tissue	Tissue Type	Invasive Ductal Carcinoma
Elisa Titer	1 - 10 µg/mL	3.58 µg/mL
Bioburden	Pass	Pass
Visible bacterial colonies must be less than 100 CFU/mL		
Source	Supernatant	Supernatant

NOTICE: This document is proprietary and its contents are the exclusive property of Merck KGaA and its affiliates. This document may not be reproduced in any form whatsoever, without prior written permission from Quality Assurance Department of the life science business of Merck.

Template Document Number: 20107867, 4.0



Document Number

Form 7.5-58-2

Title

Certificate of Analysis - Predilute

10/20/20

Quality Assurance
Rocklin, California
US

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packaging slip for additional terms and conditions of sale.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact misbranding@sial.com

NOTICE: This document is proprietary and its contents are the exclusive property of Merck KGaA and its affiliates. This document may not be reproduced in any form whatsoever, without prior written permission from Quality Assurance Department of the life science business of Merck.

Template Document Number: 20107867, 4.0



Document Number
Form 7.5-58-2

Title
Certificate of Analysis - Predilute

REVISION SHEET

Date	Version	Nature of Changes	Changed By
12/5/18	1.0	Initial Release	Chris Kaneshiro
1/24/19	2.0	Add IHC background staining pass criteria	Chris Kaneshiro

NOTICE: This document is proprietary and its contents are the exclusive property of Merck KGaA and its affiliates. This document may not be reproduced in any form whatsoever, without prior written permission from Quality Assurance Department of the life science business of Merck.

Template Document Number: 20107867, 4.0