

Certificate of Analysis

Product Name: Beta-Catenin (14) Mouse Monoclonal Antibody
Catalog Number(s): 760-4242
Roche Part Number: 05269016001
Lot Number: V0000972
Date of Manufacture: 2017-08-22
Expiration Date: 2020-07-31

| Test | Acceptance Criteria | Results |
|---|--|--|
| Functional testing on the following tissue(s): Breast | Intensity difference between sample and a reference must be less than or equal to 0.5 | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| | Background difference between sample and a reference must be less than or equal to 0.5 | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| Appearance | Transparent liquid and free of visible foreign material. | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| Bioburden | Visible bacterial colonies must be less than 100CFU/mL. | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| ELISA | 1.0-5.0 µg/mL | 1.68 µg/mL |
| Source | N/A | Supernatant |

Additional information: N/A


CELL MARQUE
 6600 Sierra College Blvd.
 Rocklin, California 95677
 (916) 746-8900 Fax: (916) 746-8989

Authorized Signature: _____

Date: _____

Quality Assurance Representative

Issued: 6/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]