



Certificate of Analysis

Product name

KIT c4800 BRAF v600 AMP/DET 24T US-IVD

Material No.

05985579190

Lot No.

F20991

Date of Manufacturing

2019-12-05

Expiry Date

2020-09

Parameter**Result****Unit**

Components

1

C48BRAF-V600MT

PASS

The visual inspection should be as follows:

3 x RXNMIX (Reaction Mix, 0.16 mL)

Natural microtube and screw cap with a natural insert,
containing a clear, colorless liquid

3 x MGAC (Magnesium Acetate, 0.15 mL)

Natural microtube and screw cap with a yellow insert,
containing a clear, colorless liquid

3 x BRAF OM (BRAF Oligo Mix, 0.13 mL)

Black microtube and screw cap with a white insert,
containing a clear liquid with a bluish/purplish tint

2 x DNA SD (DNA Specimen Diluent, 1.0 mL)

Natural microtube and screw cap with a purple insert,
containing a clear, colorless liquid

2 x BRAF MUT (BRAF Mutant Control, 0.13 mL)

Natural microtube and screw cap with a red insert,
containing a clear, colorless liquid

Issue date 2019-12-17

Page 1 of 3

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771



Certificate of Analysis

Product name

KIT c4800 BRAF v600 AMP/DET 24T US-IVD

Material No.

05985579190

Lot No.

F20991

Date of Manufacturing

2019-12-05

Expiry Date

2020-09

Parameter	Result	Unit
2 x BRAF WT (BRAF Wild Type Control, 0.13 mL) Natural microtube and screw cap with a blue insert, containing a clear, colorless liquid		
MUT Ctrl-Ch0-Ct	26.6	#InRng
MUT Ctrl-Ch1-Ct	25.8	#InRng
WT Ctrl-Ch0-Ct	1	#InRng
WT Ctrl-Ch1-Ct	34.2	#InRng
DMUT Ctrl-Mut Det	20	#InRng
DMUT Ctrl-Ch0 Ct	20	#InRng
DMUT Ctrl-Ch1 Ct	20	#InRng

This lot meets our specification

Disposition by

Issue date 2019-12-17

Page 2 of 3

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771



Certificate of Analysis

Product name

KIT c4800 BRAF v600 AMP/DET 24T US-IVD

Material No.

05985579190

Lot No.

F20991

Date of Manufacturing

2019-12-05

Expiry Date

2020-09

Craig Cawley

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2019-12-17

Page 3 of 3

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771