

CERTIFICATE OF ANALYSIS


Product Name: anti-ERG (EPR3864) Rabbit Monoclonal Primary Antibody
 Ventana REF: 790-4576
 Roche Number: 06478450001
 Batch Number: D07927

Test	Acceptance Criteria	Result
Functional Staining on qualified prostate carcinoma	Test slide staining must score 3.0 or greater on prostate tissue with a background intensity within 0.5 of the negative control slide on a scale of 0-4+	PASS
	The specific stain intensity must be within 0.5 pt. of the reference batch	PASS
	The intensity of the negative reagent control slide must be 0.5 or less	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of intrinsic particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiry date: 2014-10-28

Manufacture date: 15-OCT-2013



Quality Department **Date**