



Certificate of Analysis

Product name

KIT COBAS 4800 PIK3CA AMP/DET 24T RUO

Material No.

06523013190

Lot No.

G31757

Date of Manufacturing

2021-02-25

Expiry Date

2022-12

Parameter	Result	Unit
Components	1	
PIK3CA-MUT24T	PASS	
2 x PIK3CA MMX-1 (PIK3CA Master Mix-1, 0.48 mL) Black microtube and screw cap with a white insert, containing a clear liquid with a bluish/purplish tint		
2 x PIK3CA MMX-2 (PIK3CA Master Mix-2, 0.48 mL) Black microtube and screw cap with a gold insert, containing a clear liquid with a bluish/purplish tint		
2 x PIK3CA MMX-3 (PIK3CA Master Mix-3, 0.48 mL) Black microtube and screw cap with a teal insert, containing a clear liquid with a bluish/purplish tint		
6 x MGAC (Magnesium Acetate, 0.20 mL) Natural microtube and screw cap with a yellow insert, containing a clear, colorless liquid		
6 x PIK3CA MC (PIK3CA Mutant Control, 0.10 mL) Natural microtube and screw cap with a red insert, containing a clear, colorless liquid		
2 x DNA SD (DNA Specimen Diluent, 3.5 mL) Natural bottle and screw cap, containing a clear,		

Issue date 2021-03-11

Page 1 of 4

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771



Certificate of Analysis

Product name

KIT COBAS 4800 PIK3CA AMP/DET 24T RUO

Material No.

06523013190

Lot No.

G31757

Date of Manufacturing

2021-02-25

Expiry Date

2022-12

Parameter	Result	Unit
colorless liquid		
MMX1-MC-E542K-F0	100	%HitRa
MMX1-MC-N345K-F1	100	%HitRa
MMX1-MC-H1047X-F3	100	%HitRa
MMX1-MC-IC-F5	100	%HitRa
MMX2-MC-E545X-F0	100	%HitRa
MMX2-MC-C420R-F1	100	%HitRa
MMX2-MC-G1049R-F3	100	%HitRa
MMX2-MC-IC-F5	100	%HitRa
MMX3-MC-Q546X-F0	100	%HitRa
MMX3-MC-R88Q-F1	100	%HitRa
MMX3-MC-M1043I-F3	100	%HitRa

Issue date 2021-03-11

Page 2 of 4

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771



Certificate of Analysis

Product name

KIT COBAS 4800 PIK3CA AMP/DET 24T RUO

Material No.

06523013190

Lot No.

G31757

Date of Manufacturing

2021-02-25

Expiry Date

2022-12

Parameter	Result	Unit
MMX3-MC-IC-F5	100	%HitRa
MMX1-SD-E542K-F0	0	%HitRa
MMX1-SD-N345K-F1	0	%HitRa
MMX1-SD-H1047X-F3	0	%HitRa
MMX1-SD-IC-F5	0	%HitRa
MMX2-SD-E545X-F0	0	%HitRa
MMX2-SD-C420R-F1	0	%HitRa
MMX2-SD-G1049R-F3	0	%HitRa
MMX2-SD-IC-F5	0	%HitRa
MMX3-SD-Q546X-F0	0	%HitRa
MMX3-SD-R88Q-F1	0	%HitRa

Issue date 2021-03-11

Page 3 of 4

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771



Certificate of Analysis

Product name

KIT COBAS 4800 PIK3CA AMP/DET 24T RUO

Material No.

06523013190

Lot No.

G31757

Date of Manufacturing

2021-02-25

Expiry Date

2022-12

Parameter**Result****Unit**

MMX3-SD-M1043I-F3

0

%HitRa

MMX3-SD-IC-F5

0

%HitRa

This lot meets our specification

Disposition by

Mujib Walli

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2021-03-11

Page 4 of 4

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771