



Certificate of Analysis

Product name
CONFIRM Anti-PR (1E2) Primary Antibody

Material No
05277990001 or 05278392001

Lot No.	Date of Manufacturing	Expiry Date
E18135	21-Aug-2018	2020-08-06

Parameter	Result	Unit
------------------	---------------	-------------

Visual Inspection	Pass	
-------------------	------	--

Reagent must be free of foreign material and substantially free of particulate matter

Reagent must be clear and free of turbidity

Functional Staining	Pass	
---------------------	------	--

Breast Carcinoma Tissue

The nuclear staining will be interpreted according to the following approximate percent tumor cell staining criteria

0-50% cell staining (negative)

51-100% cell staining (positive)

The specific stain intensity must also be within 0.5 pt of the reference batch for a passing result (Scale 0-4+)

The background intensity must be within 0.5 of the negative control slide (Scale 0-4+)

The intensity of the negative reagent control slide must be 0.5 or less (Scale 0-4+)

Bioburden	Pass	
-----------	------	--

Visible bacterial colonies must be 100 or less CFU/mL

This lot meets our specification

Disposition by

JOSHUA BURGETTE

Issue date 10/Sep/2018