

Customer Address:

Customer Number : 2049990013
Customer Mat # : 06523927001
Roche Diagnostics Corporation
Accounts Payable
Po Box 50457
INDIANAPOLIS IN 46250-0417
UNITED STATES

Certificate of Analysis

Product Name : Anti-PAX8 (MRQ-50) Mouse Monoclonal Primary Antibody
Product Number : 760-4618
Batch Number : V0003402
Source Batch : 0000177868
Storage Temperature : 2°C - 8°C
Expiry Date : 2025-05-31
Quality Release Date : 2022-08-03
Date of Manufacture : 2022-06-07

| Test | Specification | Result |
|---|----------------------------------|----------------------------------|
| Appearance (Turbidity) | Free of Visible Foreign Material | Free of Visible Foreign Material |
| Appearance (Color) | Green Transparent | Green Transparent |
| Appearance (Form) | Liquid | Liquid |
| IHC Staining | Pass | Pass |
| Intensity difference between sample and a reference must be less than or equal to 0.5. | | |
| IHC Staining | Pass | Pass |
| Background difference between sample and a reference must be less than or equal to 0.5. | | |
| IHC - Tissue | Tissue Type | Ovarian Carcinoma |
| ELISA Titer | 5 - 15 µg/ml | 5 µg/ml |
| Bioburden | Pass | Pass |
| Visible bacterial colonies must be less than 100 CFU/mL. | | |
| Source | Supernatant | Supernatant |

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.



Product 760-4618
Batch V0003402



Oliver Ghattas, Quality, QC Lab Technician II

Quality Assurance

Rocklin, California

US

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.

