

### Certificate of Analysis

**Product Name:** anti-p57/Kip2 (Kp10) Mouse Monoclonal Primary Antibody  
**Catalog Number(s):** 760-4617  
**Roche Part Number:** 06523897001  
**Lot Number:** V0001066  
**Date of Manufacture:** 2017-12-12  
**Expiration Date:** 2020-11-30

Test	Acceptance Criteria	Results
Functional testing on the following tissue(s): Placenta	Intensity difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
	Background difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Appearance	Transparent liquid and free of visible foreign material.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Bioburden	Visible bacterial colonies must be less than 100CFU/mL.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
ELISA	1.0-10.0 µg/mL	1.11 µg/mL
Source	N/A	Supernatant

Additional information: N/A



Authorized Signature:  Date: 01/23/2018  
 Quality Assurance Representative

Issued: 6/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]