

Customer Address:

Customer Number : 2049990013
 Customer Mat # : 05463530001
 Roche Diagnostics Corporation
 Accounts Payable
 Po Box 50457
 INDIANAPOLIS IN 46250-0417
 UNITED STATES

Certificate of Analysis

Product Name : GCDFP-15 (EP1582Y) Rabbit Monoclonal Antibody
Product Number : 760-4386
Batch Number : V0003244
 Source Batch : 0000148239
 Storage Temperature : 2°C - 8°C
 Expiry Date : 2024-12-31
 Quality Release Date : 2022-04-08
 Date of Manufacture : 2022-01-25

| Test | Specification | Result |
|---|----------------------------------|----------------------------------|
| Appearance (Turbidity) | Free of Visible Foreign Material | Free of Visible Foreign Material |
| Appearance (Color) | Green Transparent | Green Transparent |
| Appearance (Form) | Liquid | Liquid |
| IHC Staining | Pass | Pass |
| Intensity difference between sample and a reference must be less than or equal to 0.5. | | |
| IHC Staining | Pass | Pass |
| Background difference between sample and a reference must be less than or equal to 0.5. | | |
| IHC - Tissue | Tissue Type | Breast |
| ELISA Titer | 0.01 - 1.00 µg/ml | 0.40 µg/ml |
| Bioburden | Pass | Pass |
| Visible bacterial colonies must be less than 100 CFU/mL. | | |
| Source | Supernatant | Supernatant |

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.



Product 760-4386
Batch V0003244



Oliver Ghattas, Quality, QC Lab Technician II
Quality Assurance
Rocklin, California
US

Cell Marque™ (a Sigma-Aldrich Company) warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. For further inquiries, please contact Technical/Customer Service. Purchase must determine the suitability of the product for its particular use. See website or packing slip for additional terms and conditions of sale.

