

# Certificate of Analysis



**Product Name**  
Gram Staining Kit

**Material No.**  
06890105001

**Lot No.**  
F22427

**Date of Manufacturing**  
02/OCT/2019

**Expiry Date**  
2020/09/09

**Parameter**

**Result**

**Unit**

Final Acceptance Testing  
on tissue containing Gram-positive  
and Gram-negative bacteria.

PASS

Gram-positive bacteria stain blue and dark purple.

Gram-negative bacteria stain pink to red.

Background will be yellow (Tartrazine) or green  
(Fast Green) and does not interfere with evaluation  
of Gram-positive and Gram-negative  
bacteria.

This lot meets our specification.

**Disposition By**

Guadalupe Villarreal

*Guadalupe Villarreal GV*

**Issue Date**

*04 NOV 2019*