



a member of the Roche group

CERTIFICATE OF ANALYSIS

Product Name: ALCIAN BLUE STAINING KIT

Ventana **REF** : 860-002

Roche Number: 05279186001

Batch Number: F00991

Test	Acceptance Criteria	Result
Slide Staining Test	Alcian blue stains weakly acidic mucosubstances bright blue against a pink to red background. Goblet cells which synthesize and secrete mucus stain blue. Negative control - Pink to red counterstained tissue.	PASS
Final Container Test	Staining comparable to the reference slide.	PASS

This lot met the specifications.

Expiry date: 2016-09-28

Manufacture date : 09-Apr-2015

Guadalupe Villarreal

Guadalupe Villarreal

10 APR 2015

Quality Department

Date