



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Cover letter for product:

Trade name : cobas SARS-CoV-2 Qualitative
Product code : 09446109190

The product is sold as a kit, and contains the following components:

- PASE
- DNA IC/QS/GSD
- EB
- MMX-R2
- MMX-R1

The following is an overview of the labeling elements of the kit:

GHS label elements

No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Other hazards

None known.

The product was evaluated per International Air Transport Association (IATA) specifications with the following outcome:

Not assigned by regulation


cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020
SECTION 1. IDENTIFICATION

Product name : PASE

Manufacturer or supplier's details

Company name of supplier : Roche Diagnostics Deutschland GmbH

Address : 116 Sandhoferstrasse
Mannheim, 68305
Deutschland

Telephone : +496217590

Telefax : +496217592890

E-mail address : info.dia-sds@roche.com

Emergency telephone

Im Notfall: : Werkschutzzentrale Roche +49(0)621-759-2203
Diagnostics GmbH

Giftnotruf: : Mainz +49(0)6131-19240

Recommended use of the chemical and restrictions on use
Recommended use : Laboratory chemicals
Refer to product literature for further details.
SECTION 2. HAZARDS IDENTIFICATION
GHS classification in accordance with the OSHA Hazard Communication Standard (29 CFR 1910.1200)

Not a hazardous substance or mixture.

GHS label elements

No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Other hazards

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Components

Chemical name	CAS-No.	Concentration (% w/w)
1,2,3-Propanetriol	56-81-5	>= 50 - < 70
Subtilisin	9014-01-1	< 0.1

Actual concentration is withheld as a trade secret

SECTION 4. FIRST AID MEASURES



cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- General advice : Do not leave the victim unattended.
- If inhaled : If unconscious, place in recovery position and seek medical advice.
If symptoms persist, call a physician.
- In case of skin contact : Wash off with soap and water.
- In case of eye contact : Remove contact lenses.
Protect unharmed eye.
If eye irritation persists, consult a specialist.
- If swallowed : Keep respiratory tract clear.
Do not give milk or alcoholic beverages.
Never give anything by mouth to an unconscious person.
If symptoms persist, call a physician.
- Most important symptoms and effects, both acute and delayed : None known.
- Notes to physician : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

- Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Unsuitable extinguishing media : High volume water jet
- Specific hazards during fire fighting : No information available.
- Hazardous combustion products : Carbon oxides
- Further information : Standard procedure for chemical fires.
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Special protective equipment for fire-fighters : Wear self-contained breathing apparatus for firefighting if necessary.

SECTION 6. ACCIDENTAL RELEASE MEASURES



cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
 Date of first issue: 02-26-2020

- Personal precautions, protective equipment and emergency procedures : Refer to protective measures listed in sections 7 and 8.
- Environmental precautions : If the product contaminates rivers and lakes or drains inform respective authorities.
- Methods and materials for containment and cleaning up : Wipe up with absorbent material (e.g. cloth, fleece).
Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

- Advice on protection against fire and explosion : Normal measures for preventive fire protection.
- Advice on safe handling : For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.
- Conditions for safe storage : Electrical installations / working materials must comply with the technological safety standards.
- Further information on storage conditions : See label, package insert or internal guidelines
- Materials to avoid : No materials to be especially mentioned.
- Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
1,2,3-Propanetriol	56-81-5	TWA (mist, respirable fraction)	5 mg/m ³	OSHA Z-1
		TWA (mist, total dust)	15 mg/m ³	OSHA Z-1
		TWA (Mist - total dust)	10 mg/m ³	OSHA P0
		TWA (Mist - respirable fraction)	5 mg/m ³	OSHA P0
Subtilisin	9014-01-1	IOEL	0.00006 mg/m ³	Roche Industrial Hygiene

SAFETY DATA SHEET



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

				Committee (RIHC)
		C	0.00006 mg/m3	OSHA P0
		C	0.00006 mg/m3 (crystalline active enzyme)	ACGIH
		ST	0.00006 mg/m3	NIOSH REL

Engineering measures : No data available

Personal protective equipment

Respiratory protection : General and local exhaust ventilation is recommended to maintain vapor exposures below recommended limits. Where concentrations are above recommended limits or are unknown, appropriate respiratory protection should be worn. Follow OSHA respirator regulations (29 CFR 1910.134) and use NIOSH/MSHA approved respirators. Protection provided by air purifying respirators against exposure to any hazardous chemical is limited. Use a positive pressure air supplied respirator if there is any potential for uncontrolled release, exposure levels are unknown, or any other circumstance where air purifying respirators may not provide adequate protection.

Hand protection

In case of contact through splashing:
Material : Nitrile rubber
Break through time : > 30 min
Glove thickness : > 0.11 mm

In case of full contact:
Material : butyl-rubber
Break through time : > 480 min
Glove thickness : > 0.4 mm

Remarks : Wear appropriate protective gloves to prevent skin contact.
Replace torn or punctured gloves promptly.

Eye protection : Safety glasses

Skin and body protection : Protective suit

Hygiene measures : General industrial hygiene practice.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Color : yellow



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- Odor : No data available

- Odor Threshold : No data available

- pH : 5.4 - 5.6
Concentration: 100 %

- Melting point/ range : < -85 °F / < -65 °C

- Boiling point/boiling range : 208 °F / 98 °C

- Flash point : No data available

- Evaporation rate : No data available

- Flammability (solid, gas) : The product is not flammable., Does not sustain combustion.

- Flammability (liquids) : Does not sustain combustion.

- Self-ignition : No data available

- Upper explosion limit / Upper flammability limit : No data available

- Lower explosion limit / Lower flammability limit : No data available

- Vapor pressure : No data available

- Relative vapor density : No data available

- Relative density : No data available

- Density : 1.16 g/cm³

- Solubility(ies)
 - Water solubility : completely miscible

 - Solubility in other solvents : No data available

- Partition coefficient: n-octanol/water : No data available

- Autoignition temperature : No data available

- Decomposition temperature : No data available

- Viscosity
 - Viscosity, dynamic : No data available

 - Viscosity, kinematic : No data available



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Explosive properties : Not explosive

Oxidizing properties : The substance or mixture is not classified as oxidizing.

Particle characteristics
Particle Size Distribution : Not applicable

SECTION 10. STABILITY AND REACTIVITY

Reactivity : No dangerous reaction known under conditions of normal use.

Chemical stability : Stable under normal conditions.

Possibility of hazardous reactions : Stable under recommended storage conditions.
No hazards to be specially mentioned.

Conditions to avoid : No data available

Incompatible materials : Not applicable

Hazardous decomposition products : No hazardous decomposition products are known.

SECTION 11. TOXICOLOGICAL INFORMATION

Acute toxicity

Not classified due to lack of data.

Components:

1,2,3-Propanetriol:

Acute oral toxicity : LC50 (Mouse): 11,500 mg/kg

Acute inhalation toxicity : LC50 (Rat, male): 275000 mg/m³
Exposure time: 7 h
Test atmosphere: vapor
GLP: no
Assessment: The component/mixture is minimally toxic after short term inhalation.

Acute dermal toxicity : LD50 (Guinea pig, male and female): 56,750 mg/kg
GLP: no

Subtilisin:

Acute oral toxicity : LD50 Oral (Rat, male and female): 1,800 mg/kg
Method: OECD Test Guideline 401
GLP: yes



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Skin corrosion/irritation

Not classified due to lack of data.

Components:

1,2,3-Propanetriol:

Species : Rabbit
Exposure time : 24 h
Result : No skin irritation
GLP : no

Subtilisin:

Species : Rabbit
Method : OECD Test Guideline 404
Result : Irritating to skin.
GLP : yes

Serious eye damage/eye irritation

Not classified due to lack of data.

Components:

1,2,3-Propanetriol:

Species : Rabbit
Result : No eye irritation
Exposure time : 7 d
GLP : no

Subtilisin:

Species : Rabbit
Result : Risk of serious damage to eyes.
Method : OECD Test Guideline 405
GLP : yes

Respiratory or skin sensitization

Skin sensitization

Not classified due to lack of data.

Respiratory sensitization

Not classified due to lack of data.

Components:

1,2,3-Propanetriol:

Assessment : Mild eye irritant, Mild respiratory irritant, No skin irritation



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Subtilisin:

Assessment : May cause sensitization by inhalation.

Germ cell mutagenicity

Not classified due to lack of data.

Components:

1,2,3-Propanetriol:

Genotoxicity in vitro : Test Type: Microbial mutagenesis assay (Ames test)
Test system: Salmonella typhimurium
Metabolic activation: with and without metabolic activation
Result: negative
GLP: No information available.

Test Type: In vitro mammalian cell gene mutation test
Test system: Chinese hamster ovary cells
Metabolic activation: with and without metabolic activation
Method: OECD Test Guideline 476
Result: negative
GLP: No information available.

Subtilisin:

Genotoxicity in vitro : Test Type: In vitro mammalian cell gene mutation test
Test system: mouse lymphoma cells
Result: negative

Test Type: Microbial mutagenesis assay (Ames test)
Test system: Escherichia coli
Result: negative

Test Type: Microbial mutagenesis assay (Ames test)
Test system: Salmonella typhimurium
Result: negative

Test Type: Chromosome aberration test in vitro
Test system: Human lymphocytes
Result: negative

Carcinogenicity

Not classified due to lack of data.



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Components:

1,2,3-Propanetriol:

Species : Rat, male and female
Application Route : Oral
Exposure time : 2 Years
GLP : No information available.
Remarks : No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

Reproductive toxicity

Not classified due to lack of data.

Components:

1,2,3-Propanetriol:

Effects on fertility : Test Type: Two-generation study
Species: Rat, male and female
Application Route: Oral
Dose: 2000 mg/kg bw/day
Fertility: NOAEL: 2,000 mg/kg body weight
GLP: no

Effects on fetal development : Species: Rabbit, female
Application Route: Oral
Dose: 11.8, 54.8, 254.5, 1180 mg/kg bw/day
Duration of Single Treatment: 29 d
Developmental Toxicity: NOAEL: 1,180 mg/kg bw/day
GLP: no

STOT-single exposure

Not classified due to lack of data.

Components:

Subtilisin:

Routes of exposure : Inhalation
Target Organs : Lungs
Assessment : May cause respiratory irritation.

STOT-repeated exposure

Not classified due to lack of data.

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020**Repeated dose toxicity****Components:****1,2,3-Propanetriol:**

Species : Rat, male and female
 NOAEL : 4580 mg/kg
 NOAEL : 4,580 mg/kg
 Application Route : Oral
 Exposure time : 90 d
 Number of exposures : daily
 Dose : 4580 - 25,800 mg/kg/day
 GLP : no

Species : Rat, male and female
 Application Route : Inhalation
 Test atmosphere : dust/mist
 Exposure time : 13 Weeks
 Number of exposures : 6 hours/day, 5 days/week
 Dose : 33, 165 and 660 mg/m³
 GLP : No information available.

Species : Rat
 NOAEL : 5040 mg/kg
 NOAEL : 5,040 mg/kg
 Application Route : Dermal
 Exposure time : 45 Weeks
 Number of exposures : 8 hours/day, 5 days/week
 Dose : 0.5-4.0 ml/kg
 GLP : no

Repeated dose toxicity - Assessment : Mild eye irritant, Mild respiratory irritant, No skin irritation

Aspiration toxicity

Not classified due to lack of data.

Further information**Product:**

Remarks : No data available

SECTION 12. ECOLOGICAL INFORMATION
Ecotoxicity**Product:****Ecotoxicology Assessment**

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020**Components:****1,2,3-Propanetriol:**

Toxicity to fish : LC50 (Oncorhynchus mykiss (rainbow trout)): 54,000 mg/l
 End point: mortality
 Exposure time: 96 h
 Test Type: static test
 GLP: no

Toxicity to daphnia and other aquatic invertebrates : LC50 (Daphnia magna (Water flea)): 1,955 mg/l
 End point: mortality
 Exposure time: 48 h
 Test Type: static test
 Analytical monitoring: no
 GLP: no

Toxicity to algae/aquatic plants : (Scenedesmus quadricauda (Green algae)): > 10,000 mg/l
 End point: Growth rate
 Exposure time: 8 d
 Test Type: static test
 GLP: no

Toxicity to microorganisms : EC50 (Pseudomonas putida): > 10,000 mg/l
 End point: Growth rate
 Exposure time: 16 h
 Test Type: static test
 GLP: No information available.

Ecotoxicology Assessment

Acute aquatic toxicity : This product has no known ecotoxicological effects.

Chronic aquatic toxicity : This product has no known ecotoxicological effects.

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Subtilisin:

Toxicity to fish : LC50 (Oncorhynchus mykiss (rainbow trout)): 8.2 mg/l
 Exposure time: 96 h
 Test Type: semi-static test
 Method: OECD Test Guideline 203
 GLP: yes

Toxicity to daphnia and other aquatic invertebrates : EC50 (Daphnia magna (Water flea)): 0.306 mg/l
 End point: Immobilization
 Exposure time: 48 h
 Test Type: semi-static test
 Method: OECD Test Guideline 202
 Remarks: average measured concentration

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- NOEC (*Daphnia magna* (Water flea)): 0.066 mg/l
End point: Immobilization
Exposure time: 48 h
Test Type: semi-static test
Method: OECD Test Guideline 202
Remarks: average measured concentration
- Toxicity to algae/aquatic plants : ErC50 (*Raphidocelis subcapitata* (freshwater green alga)): 0.830 mg/l
End point: Growth rate
Exposure time: 72 h
Test Type: static test
Method: OECD Test Guideline 201
GLP: yes
Remarks: measured initial concentration
- NOEC (*Raphidocelis subcapitata* (freshwater green alga)): 0.566 mg/l
End point: Growth rate
Exposure time: 72 h
Test Type: static test
Method: OECD Test Guideline 201
GLP: yes
Remarks: measured initial concentration
- Toxicity to fish (Chronic toxicity) : NOEC (*Pimephales promelas* (fathead minnow)): 0.042 mg/l
End point: mortality
Exposure time: 32 d
Test Type: Fish early-life stage (FELS) toxicity test (OECD 210)
Method: OECD Test Guideline 210
GLP: yes
Remarks: average measured concentration

Persistence and degradability
Components:
1,2,3-Propanetriol:

- Biodegradability : aerobic
Inoculum: activated sludge
Concentration: 226 mg/l
Result: Readily biodegradable.
Biodegradation: 94 %
Exposure time: 24 h
GLP: no

Subtilisin:

- Biodegradability : aerobic
Result: Readily biodegradable.
Biodegradation: 100 %
Exposure time: 29 d
Method: OECD Test Guideline 301B



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

GLP: yes

Bioaccumulative potential

Components:

1,2,3-Propanetriol:

Partition coefficient: n-octanol/water : log Pow: -1.75 (77 °F / 25 °C)
pH: 7.4
Method: OECD Test Guideline 107
GLP: no

Subtilisin:

Partition coefficient: n-octanol/water : log Pow: -3.1 (77 °F / 25 °C)
pH: 9.2
Method: OECD Test Guideline 107
GLP: no

Mobility in soil

No data available

Other adverse effects

Product:

Ozone-Depletion Potential : Regulation: 40 CFR Protection of Environment; Part 82
Protection of Stratospheric Ozone - CAA Section 602 Class I
Substances
Remarks: This product neither contains, nor was
manufactured with a Class I or Class II ODS as defined by the
U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A +
B).

Additional ecological information : No data available

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : Can be disposed as waste water, when in compliance with
local regulations.

Contaminated packaging : Empty containers should be taken to an approved waste
handling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR



cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
 Date of first issue: 02-26-2020

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

Domestic regulation

49 CFR

Not regulated as a dangerous good

Special precautions for user

Remarks : Not dangerous goods in the meaning of ADR/RID, ADN, IMDG-Code, ICAO/IATA-DGR

SECTION 15. REGULATORY INFORMATION

CERCLA Reportable Quantity

This material does not contain any components with a CERCLA RQ.

SARA 304 Extremely Hazardous Substances Reportable Quantity

This material does not contain any components with a section 304 EHS RQ.

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

This material does not contain any components with a section 302 EHS TPQ.

SARA 311/312 Hazards : No SARA Hazards

SARA 313 : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Clean Air Act

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

The following chemical(s) are listed under the U.S. Clean Air Act Section 111 SOCM I

Intermediate or Final VOC's (40 CFR 60.489):

1,2,3-Propanetriol	56-81-5	>= 50 - < 70 %
--------------------	---------	----------------

Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

This product does not contain any priority pollutants related to the U.S. Clean Water Act



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

US State Regulations

Massachusetts Right To Know

1,2,3-Propanetriol 56-81-5

Pennsylvania Right To Know

1,2,3-Propanetriol 56-81-5
Water 7732-18-5

Maine Chemicals of High Concern

Product does not contain any listed chemicals

Vermont Chemicals of High Concern

Product does not contain any listed chemicals

Washington Chemicals of High Concern

Product does not contain any listed chemicals

California Permissible Exposure Limits for Chemical Contaminants

1,2,3-Propanetriol 56-81-5

The ingredients of this product are reported in the following inventories:

- AIIC : All components are listed on the inventory, regulatory obligations/restrictions apply
- DSL : This product contains the following components that are not on the Canadian DSL nor NDSL.

Tetrasodium ethylenediaminetetraacetate dihydrate
- NZIoC : On the inventory, or in compliance with the inventory
- ENCS : Not in compliance with the inventory
- ISHL : On the inventory, or in compliance with the inventory
- KECI : On the inventory, or in compliance with the inventory
- PICCS : On the inventory, or in compliance with the inventory
- IECSC : On the inventory, or in compliance with the inventory
- TCSI : On the inventory, or in compliance with the inventory
- TSCA : Product contains substance(s) not listed on TSCA inventory.
- TECI : Not in compliance with the inventory

TSCA list

No substances are subject to a Significant New Use Rule.

No substances are subject to TSCA 12(b) export notification requirements.

cobas SARS-CoV-2 Qualitative

Version
4.0


Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 16. OTHER INFORMATION

Further information

NFPA 704:



HMIS® IV:

HEALTH	/	0
FLAMMABILITY		0
PHYSICAL HAZARD		0

HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. The "*" represents a chronic hazard, while the "/" represents the absence of a chronic hazard.

Full text of other abbreviations

- ACGIH : USA. ACGIH Threshold Limit Values (TLV)
- NIOSH REL : USA. NIOSH Recommended Exposure Limits
- OSHA P0 : USA. Table Z-1-A Limits for Air Contaminants (1989 vacated values)
- OSHA Z-1 : USA. Occupational Exposure Limits (OSHA) - Table Z-1 Limits for Air Contaminants
- ACGIH / C : Ceiling limit
- NIOSH REL / ST : STEL - 15-minute TWA exposure that should not be exceeded at any time during a workday
- OSHA P0 / TWA : 8-hour time weighted average
- OSHA P0 / C : Ceiling limit
- OSHA Z-1 / TWA : 8-hour time weighted average

AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO -



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 08-20-2025

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

US / Z8 / 2404



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- If inhaled : If unconscious, place in recovery position and seek medical advice.
If symptoms persist, call a physician.
- In case of skin contact : Wash off with soap and water.
- In case of eye contact : Remove contact lenses.
Protect unharmed eye.
If eye irritation persists, consult a specialist.
- If swallowed : Keep respiratory tract clear.
Do not give milk or alcoholic beverages.
Never give anything by mouth to an unconscious person.
If symptoms persist, call a physician.
- Most important symptoms and effects, both acute and delayed : None known.
- Notes to physician : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

- Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Unsuitable extinguishing media : High volume water jet
- Specific hazards during fire fighting : No information available.
- Hazardous combustion products : No hazardous combustion products are known
- Further information : Standard procedure for chemical fires.
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Special protective equipment for fire-fighters : Wear self-contained breathing apparatus for firefighting if necessary.

SECTION 6. ACCIDENTAL RELEASE MEASURES

- Personal precautions, protective equipment and : Refer to protective measures listed in sections 7 and 8.



cobas SARS-CoV-2 Qualitative

Version 4.0	Revision Date: 08-20-2025	Date of last issue: 08-20-2025 Date of first issue: 02-26-2020
----------------	------------------------------	---

emergency procedures

Environmental precautions : No special environmental precautions required.

Methods and materials for containment and cleaning up : Wipe up with absorbent material (e.g. cloth, fleece).
Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

Advice on protection against fire and explosion : Normal measures for preventive fire protection.

Advice on safe handling : For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.

Conditions for safe storage : Electrical installations / working materials must comply with the technological safety standards.

Further information on storage conditions : See label, package insert or internal guidelines

Materials to avoid : No materials to be especially mentioned.

Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Contains no substances with occupational exposure limit values.

Engineering measures : No data available

Personal protective equipment

Hand protection

		In case of contact through splashing:
Material	:	Nitrile rubber
Break through time	:	> 30 min
Glove thickness	:	> 0.11 mm

		In case of full contact:
Material	:	butyl-rubber
Break through time	:	> 480 min
Glove thickness	:	> 0.4 mm

Remarks : Wear appropriate protective gloves to prevent skin contact.
Replace torn or punctured gloves promptly.

Eye protection : Safety glasses



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Skin and body protection : Protective suit
Hygiene measures : General industrial hygiene practice.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Color : No data available

Odor : No data available

Odor Threshold : No data available

pH : 8.0

Melting point/ range : No data available

Boiling point/boiling range : No data available

Flash point : does not flash

Evaporation rate : No data available

Flammability (solid, gas) : The product is not flammable., Does not sustain combustion.

Flammability (liquids) : Does not sustain combustion.

Self-ignition : No data available

Upper explosion limit / Upper flammability limit : No data available

Lower explosion limit / Lower flammability limit : No data available

Vapor pressure : No data available

Relative vapor density : No data available

Relative density : No data available

Density : 1.006 g/cm³ (68 °F / 20 °C, 990 - 1,030 hPa)

Solubility(ies)
Water solubility : completely miscible



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Solubility in other solvents	:	No data available
Partition coefficient: n-octanol/water	:	No data available
Autoignition temperature	:	No data available
Decomposition temperature	:	No data available
Viscosity		
Viscosity, dynamic	:	No data available
Viscosity, kinematic	:	No data available
Explosive properties	:	Not explosive
Oxidizing properties	:	The substance or mixture is not classified as oxidizing.
Particle characteristics		
Particle Size Distribution	:	Not applicable

SECTION 10. STABILITY AND REACTIVITY

Reactivity	:	No dangerous reaction known under conditions of normal use. No dangerous reaction known under conditions of normal use.
Chemical stability	:	Stable under normal conditions.
Possibility of hazardous reactions	:	Stable under recommended storage conditions. No hazards to be specially mentioned.
Conditions to avoid	:	No data available
Incompatible materials	:	No data available Not applicable
Hazardous decomposition products	:	No hazardous decomposition products are known.

SECTION 11. TOXICOLOGICAL INFORMATION

Acute toxicity

Not classified due to lack of data.

Skin corrosion/irritation

Not classified due to lack of data.

Serious eye damage/eye irritation

Not classified due to lack of data.



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Respiratory or skin sensitization

Skin sensitization

Not classified due to lack of data.

Respiratory sensitization

Not classified due to lack of data.

Germ cell mutagenicity

Not classified due to lack of data.

Carcinogenicity

Not classified due to lack of data.

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

Reproductive toxicity

Not classified due to lack of data.

STOT-single exposure

Not classified due to lack of data.

STOT-repeated exposure

Not classified due to lack of data.

Aspiration toxicity

Not classified due to lack of data.

Further information

Product:

Remarks : No data available

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Product:

Ecotoxicology Assessment

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Persistence and degradability

No data available



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

Product:

Ozone-Depletion Potential : Regulation: 40 CFR Protection of Environment; Part 82
Protection of Stratospheric Ozone - CAA Section 602 Class I
Substances
Remarks: This product neither contains, nor was
manufactured with a Class I or Class II ODS as defined by the
U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A +
B).

Additional ecological : No data available
information

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : Can be disposed as waste water, when in compliance with
local regulations.

Contaminated packaging : Empty containers should be taken to an approved waste
handling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

Domestic regulation

49 CFR

Not regulated as a dangerous good

Special precautions for user

Remarks : Not dangerous goods in the meaning of ADR/RID, ADN,
IMDG-Code, ICAO/IATA-DGR



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 15. REGULATORY INFORMATION

CERCLA Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 304 Extremely Hazardous Substances Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

This material does not contain any components with a section 302 EHS TPQ.

SARA 311/312 Hazards : No SARA Hazards

SARA 313 : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Clean Air Act

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 111 SOCM Intermediate or Final VOC's (40 CFR 60.489).

Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

This product does not contain any priority pollutants related to the U.S. Clean Water Act

US State Regulations

Massachusetts Right To Know

Sodium azide (Na(N3)) 26628-22-8

Pennsylvania Right To Know

Water 7732-18-5
Sodium azide (Na(N3)) 26628-22-8

Maine Chemicals of High Concern

Product does not contain any listed chemicals

Vermont Chemicals of High Concern

Product does not contain any listed chemicals

Washington Chemicals of High Concern

Product does not contain any listed chemicals

The ingredients of this product are reported in the following inventories:



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- AIIC : On the inventory, or in compliance with the inventory
- DSL : All components of this product are on the Canadian DSL
- NZIoC : On the inventory, or in compliance with the inventory
- ENCS : Not in compliance with the inventory
- ISHL : Not in compliance with the inventory
- KECI : Not in compliance with the inventory
- PICCS : On the inventory, or in compliance with the inventory
- IECSC : On the inventory, or in compliance with the inventory
- TCSI : On the inventory, or in compliance with the inventory
- TSCA : All substances listed as active on the TSCA inventory
- TECI : Not in compliance with the inventory

TSCA list

No substances are subject to a Significant New Use Rule.

No substances are subject to TSCA 12(b) export notification requirements.

SECTION 16. OTHER INFORMATION

Further information


cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

NFPA 704:



HMIS® IV:

HEALTH	/	0
FLAMMABILITY		0
PHYSICAL HAZARD		0

HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. The "*" represents a chronic hazard, while the "/" represents the absence of a chronic hazard.

Full text of other abbreviations

AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECl - Thailand

SAFETY DATA SHEET



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 08-20-2025

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

US / Z8 / 2404



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 1. IDENTIFICATION

Product name : EB

Manufacturer or supplier's details

Company name of supplier : Roche Diagnostics Deutschland GmbH

Address : 116 Sandhoferstrasse
Mannheim, 68305
Deutschland

Telephone : +496217590

Telefax : +496217592890

E-mail address : info.dia-sds@roche.com

Emergency telephone

Im Notfall: : Werkschutzzentrale Roche +49(0)621-759-2203
Diagnostics GmbH

Giftnotruf: : Mainz +49(0)6131-19240

Recommended use of the chemical and restrictions on use

Recommended use : Laboratory chemicals
Refer to product literature for further details.

SECTION 2. HAZARDS IDENTIFICATION

GHS classification in accordance with the OSHA Hazard Communication Standard (29 CFR 1910.1200)

Not a hazardous substance or mixture.

GHS label elements

No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Other hazards

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Components

No hazardous ingredients

SECTION 4. FIRST AID MEASURES

General advice : Do not leave the victim unattended.



cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- If inhaled : If unconscious, place in recovery position and seek medical advice.
If symptoms persist, call a physician.
- In case of skin contact : Wash off with soap and water.
- In case of eye contact : Remove contact lenses.
Protect unharmed eye.
If eye irritation persists, consult a specialist.
- If swallowed : Keep respiratory tract clear.
Do not give milk or alcoholic beverages.
Never give anything by mouth to an unconscious person.
If symptoms persist, call a physician.
- Most important symptoms and effects, both acute and delayed : None known.
- Notes to physician : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

- Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Unsuitable extinguishing media : High volume water jet
- Specific hazards during fire fighting : No information available.
- Hazardous combustion products : No hazardous combustion products are known
- Further information : Standard procedure for chemical fires.
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Special protective equipment for fire-fighters : Wear self-contained breathing apparatus for firefighting if necessary.

SECTION 6. ACCIDENTAL RELEASE MEASURES

- Personal precautions, protective equipment and : Refer to protective measures listed in sections 7 and 8.



cobas SARS-CoV-2 Qualitative

Version 4.0	Revision Date: 08-20-2025	Date of last issue: 08-20-2025 Date of first issue: 02-26-2020
----------------	------------------------------	---

emergency procedures

Environmental precautions : No special environmental precautions required.

Methods and materials for containment and cleaning up : Wipe up with absorbent material (e.g. cloth, fleece).
Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

Advice on protection against fire and explosion : Normal measures for preventive fire protection.

Advice on safe handling : For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.

Conditions for safe storage : Electrical installations / working materials must comply with the technological safety standards.

Further information on storage conditions : See label, package insert or internal guidelines

Materials to avoid : No materials to be especially mentioned.

Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Contains no substances with occupational exposure limit values.

Engineering measures : No data available

Personal protective equipment

Respiratory protection : No personal respiratory protective equipment normally required.

Hand protection

In case of contact through splashing:
 Material : Nitrile rubber
 Break through time : > 30 min
 Glove thickness : > 0.11 mm

In case of full contact:
 Material : butyl-rubber
 Break through time : > 480 min
 Glove thickness : > 0.4 mm



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Remarks : Wear appropriate protective gloves to prevent skin contact.
Replace torn or punctured gloves promptly.

Eye protection : Safety glasses

Skin and body protection : Protective suit

Hygiene measures : General industrial hygiene practice.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Color : colorless

Odor : No data available

Odor Threshold : No data available

pH : 8.3 - 8.7
Concentration: 100 %

Melting point/ range : 32.97 °F / 0.54 °C

Boiling point/boiling range : 223 °F / 106 °C

Flash point : does not flash

Evaporation rate : No data available

Flammability (solid, gas) : The product is not flammable., Does not sustain combustion.

Flammability (liquids) : Does not sustain combustion.

Self-ignition : No data available

Upper explosion limit / Upper flammability limit : No data available

Lower explosion limit / Lower flammability limit : No data available

Vapor pressure : No data available

Relative vapor density : No data available

Relative density : No data available

Density : 1 g/cm³



cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- Solubility(ies)
 - Water solubility : completely miscible
 - Solubility in other solvents : No data available
- Partition coefficient: n-octanol/water : No data available
- Autoignition temperature : No data available
- Decomposition temperature : No data available
- Viscosity
 - Viscosity, dynamic : No data available
 - Viscosity, kinematic : No data available
- Explosive properties : Not explosive
- Oxidizing properties : The substance or mixture is not classified as oxidizing.
- Particle characteristics
 - Particle Size Distribution : Not applicable

SECTION 10. STABILITY AND REACTIVITY

- Reactivity : No dangerous reaction known under conditions of normal use.
- Chemical stability : Stable under normal conditions.
- Possibility of hazardous reactions : Stable under recommended storage conditions.
No hazards to be specially mentioned.
- Conditions to avoid : No data available
- Incompatible materials : Not applicable
- Hazardous decomposition products : No hazardous decomposition products are known.

SECTION 11. TOXICOLOGICAL INFORMATION

- Acute toxicity**
Not classified due to lack of data.
- Skin corrosion/irritation**
Not classified due to lack of data.
- Serious eye damage/eye irritation**
Not classified due to lack of data.



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Respiratory or skin sensitization

Skin sensitization

Not classified due to lack of data.

Respiratory sensitization

Not classified due to lack of data.

Germ cell mutagenicity

Not classified due to lack of data.

Carcinogenicity

Not classified due to lack of data.

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

Reproductive toxicity

Not classified due to lack of data.

STOT-single exposure

Not classified due to lack of data.

STOT-repeated exposure

Not classified due to lack of data.

Aspiration toxicity

Not classified due to lack of data.

Further information

Product:

Remarks : No data available

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Product:

Ecotoxicology Assessment

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Persistence and degradability

No data available



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

Product:

Ozone-Depletion Potential : Regulation: 40 CFR Protection of Environment; Part 82
Protection of Stratospheric Ozone - CAA Section 602 Class I
Substances
Remarks: This product neither contains, nor was
manufactured with a Class I or Class II ODS as defined by the
U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A +
B).

Additional ecological : No data available
information

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : Can be disposed as waste water, when in compliance with
local regulations.

Contaminated packaging : Empty containers should be taken to an approved waste
handling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

Domestic regulation

49 CFR

Not regulated as a dangerous good

Special precautions for user

Remarks : Not dangerous goods in the meaning of ADR/RID, ADN,
IMDG-Code, ICAO/IATA-DGR



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 15. REGULATORY INFORMATION

CERCLA Reportable Quantity

This material does not contain any components with a CERCLA RQ.

SARA 304 Extremely Hazardous Substances Reportable Quantity

This material does not contain any components with a section 304 EHS RQ.

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

This material does not contain any components with a section 302 EHS TPQ.

SARA 311/312 Hazards : No SARA Hazards

SARA 313 : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Clean Air Act

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 111 SOCM Intermediate or Final VOC's (40 CFR 60.489).

Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

This product does not contain any priority pollutants related to the U.S. Clean Water Act

US State Regulations

Massachusetts Right To Know

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know

Water 7732-18-5

Maine Chemicals of High Concern

Benzoic acid, 4-hydroxy-, methyl ester 99-76-3

Vermont Chemicals of High Concern

Benzoic acid, 4-hydroxy-, methyl ester 99-76-3

Washington Chemicals of High Concern

Benzoic acid, 4-hydroxy-, methyl ester 99-76-3

The ingredients of this product are reported in the following inventories:

AIIC : On the inventory, or in compliance with the inventory



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- DSL : All components of this product are on the Canadian DSL
- NZIoC : On the inventory, or in compliance with the inventory
- ENCS : On the inventory, or in compliance with the inventory
- ISHL : On the inventory, or in compliance with the inventory
- KECI : On the inventory, or in compliance with the inventory
- PICCS : On the inventory, or in compliance with the inventory
- IECSC : On the inventory, or in compliance with the inventory
- TCSI : On the inventory, or in compliance with the inventory
- TSCA : All substances listed as active on the TSCA inventory
- TECI : On the inventory, or in compliance with the inventory

TSCA list

No substances are subject to a Significant New Use Rule.

No substances are subject to TSCA 12(b) export notification requirements.

SECTION 16. OTHER INFORMATION

Further information


cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

NFPA 704:



HMIS® IV:

HEALTH	/	0
FLAMMABILITY		0
PHYSICAL HAZARD		0

HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. The "*" represents a chronic hazard, while the "/" represents the absence of a chronic hazard.

Full text of other abbreviations

AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECl - Thailand



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 08-20-2025

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

US / Z8 / 2404

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020
SECTION 1. IDENTIFICATION

Product name : MMX-R2

Manufacturer or supplier's details

Company name of supplier : Roche Diagnostics Deutschland GmbH

Address : 116 Sandhoferstrasse
Mannheim, 68305
Deutschland

Telephone : +496217590

Telefax : +496217592890

E-mail address : info.dia-sds@roche.com

Emergency telephone

Im Notfall: : Werkschutzzentrale Roche +49(0)621-759-2203
Diagnostics GmbH

Giftnotruf: : Mainz +49(0)6131-19240

Recommended use of the chemical and restrictions on use
Recommended use : Laboratory chemicals
Refer to product literature for further details.
SECTION 2. HAZARDS IDENTIFICATION
GHS classification in accordance with the OSHA Hazard Communication Standard (29 CFR 1910.1200)

Not a hazardous substance or mixture.

GHS label elements

No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Other hazards

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Components

Chemical name	CAS-No.	Concentration (% w/w)
Methane, 1,1'-sulfinylbis-	67-68-5	>= 10 - < 20
1,2,3-Propanetriol	56-81-5	>= 10 - < 20
Acetic acid, potassium salt (1:1)	127-08-2	>= 1 - < 5
DNA-dependent DNA polymerase	9012-90-2	< 0.1

Actual concentration is withheld as a trade secret



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 4. FIRST AID MEASURES

- General advice : Do not leave the victim unattended.
- If inhaled : If unconscious, place in recovery position and seek medical advice.
If symptoms persist, call a physician.
- In case of skin contact : Wash off with soap and water.
- In case of eye contact : Remove contact lenses.
Protect unharmed eye.
If eye irritation persists, consult a specialist.
- If swallowed : Keep respiratory tract clear.
Do not give milk or alcoholic beverages.
Never give anything by mouth to an unconscious person.
If symptoms persist, call a physician.
- Most important symptoms and effects, both acute and delayed : None known.
- Notes to physician : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

- Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Unsuitable extinguishing media : High volume water jet
- Specific hazards during fire fighting : No information available.
- Hazardous combustion products : Carbon oxides
Sulfur oxides
- Further information : Standard procedure for chemical fires.
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Special protective equipment for fire-fighters : Wear self-contained breathing apparatus for firefighting if necessary.



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 6. ACCIDENTAL RELEASE MEASURES

- Personal precautions, protective equipment and emergency procedures : Refer to protective measures listed in sections 7 and 8.
- Environmental precautions : If the product contaminates rivers and lakes or drains inform respective authorities.
- Methods and materials for containment and cleaning up : Wipe up with absorbent material (e.g. cloth, fleece).
Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

- Advice on protection against fire and explosion : Normal measures for preventive fire protection.
- Advice on safe handling : For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.
- Conditions for safe storage : Electrical installations / working materials must comply with the technological safety standards.
- Further information on storage conditions : See label, package insert or internal guidelines
- Materials to avoid : No materials to be especially mentioned.
- Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
Methane, 1,1'-sulfinylbis-	67-68-5	TWA	250 ppm	US WEEL
1,2,3-Propanetriol	56-81-5	TWA (mist, respirable fraction)	5 mg/m3	OSHA Z-1
		TWA (mist, total dust)	15 mg/m3	OSHA Z-1
		TWA (Mist - total dust)	10 mg/m3	OSHA P0
		TWA (Mist - respirable fraction)	5 mg/m3	OSHA P0

SAFETY DATA SHEET



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

DNA-dependent DNA polymerase	9012-90-2	IOEL	0.00006 mg/m3	Roche Industrial Hygiene Committee (RIHC)
------------------------------	-----------	------	---------------	---

Engineering measures : No data available

Personal protective equipment

Hand protection

In case of full contact:
Material : Nitrile rubber
Break through time : > 480 min
Glove thickness : > 0.4 mm

In case of contact through splashing:
Material : Nitrile rubber
Break through time : 60 min
Glove thickness : 0.11 - 0.2 mm

Remarks : Wear appropriate protective gloves to prevent skin contact. Replace torn or punctured gloves promptly.
Eye protection : Safety glasses
Skin and body protection : Protective suit
Hygiene measures : General industrial hygiene practice.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid
Color : No data available
Odor : No data available
Odor Threshold : No data available
pH : 8.1
Melting point/ range : No data available
Boiling point/boiling range : No data available
Flash point : does not flash



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- Evaporation rate : No data available
- Flammability (solid, gas) : Does not sustain combustion.
- Flammability (liquids) : Does not sustain combustion.
- Self-ignition : Not applicable
- Upper explosion limit / Upper flammability limit : No data available
- Lower explosion limit / Lower flammability limit : No data available
- Vapor pressure : No data available
- Relative vapor density : No data available
- Relative density : No data available
- Density : 1.07 g/cm³
- Solubility(ies)
 - Water solubility : completely miscible
 - Solubility in other solvents : No data available
- Partition coefficient: n-octanol/water : No data available
- Autoignition temperature : No data available
- Decomposition temperature : No data available
- Viscosity
 - Viscosity, dynamic : No data available
 - Viscosity, kinematic : No data available
- Explosive properties : Not explosive
- Oxidizing properties : The substance or mixture is not classified as oxidizing.
- Particle characteristics
 - Particle Size Distribution : Not applicable

SECTION 10. STABILITY AND REACTIVITY

- Reactivity : No dangerous reaction known under conditions of normal use.
- Chemical stability : Stable under normal conditions.

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Possibility of hazardous reactions	:	Stable under recommended storage conditions. No hazards to be specially mentioned.
Conditions to avoid	:	No data available
Incompatible materials	:	Strong oxidizing agents Not applicable
Hazardous decomposition products	:	No data available No hazardous decomposition products are known.

SECTION 11. TOXICOLOGICAL INFORMATION
Acute toxicity

Not classified due to lack of data.

Product:Acute oral toxicity : Acute toxicity estimate: > 5,000 mg/kg
Method: Calculation method**Components:****Methane, 1,1'-sulfinylbis-:**Acute oral toxicity : LD50 (Rat, male and female): 28,300 mg/kg
Method: OECD Test Guideline 401
GLP: noAcute inhalation toxicity : LC0 (Rat, male and female): > 5.33 mg/l
Exposure time: 4 h
Test atmosphere: vapor
Method: OECD Test Guideline 403
GLP: yesAcute dermal toxicity : LD50 Dermal (Rat, male and female): 40,000 mg/kg
GLP: no**1,2,3-Propanetriol:**

Acute oral toxicity : LC50 (Mouse): 11,500 mg/kg

Acute inhalation toxicity : LC50 (Rat, male): 275000 mg/m³
Exposure time: 7 h
Test atmosphere: vapor
GLP: no
Assessment: The component/mixture is minimally toxic after short term inhalation.

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Acute dermal toxicity : LD50 (Guinea pig, male and female): 56,750 mg/kg
GLP: no

Acetic acid, potassium salt (1:1):

Acute oral toxicity : LD50 Oral (Rat): 3,250 mg/kg

Acute inhalation toxicity : Acute toxicity estimate: > 30 mg/l
Test atmosphere: dust/mist
Method: Expert judgment

LC50 (Rat): > 5.6 mg/l
Test atmosphere: dust/mist
Method: OECD Test Guideline 403
Assessment: The substance or mixture has no acute inhalation toxicity
Remarks: No mortality observed at this dose.

Acute dermal toxicity : Acute toxicity estimate: > 5,000 mg/kg
Method: Expert judgment

Skin corrosion/irritation

Not classified due to lack of data.

Components:
Methane, 1,1'-sulfinylbis-:

Species : Rabbit
Exposure time : 4 h
Method : OECD Test Guideline 404
GLP : yes
Remarks : Mild skin irritation

1,2,3-Propanetriol:

Species : Rabbit
Exposure time : 24 h
Result : No skin irritation
GLP : no

Acetic acid, potassium salt (1:1):

Species : Rabbit
Method : OECD Test Guideline 404
Result : No skin irritation

cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Serious eye damage/eye irritation

Not classified due to lack of data.

Components:
Methane, 1,1'-sulfinylbis-:

Species : Rabbit
Exposure time : 24 h
Method : OECD Test Guideline 405
GLP : No information available.
Remarks : Mild eye irritation

1,2,3-Propanetriol:

Species : Rabbit
Result : No eye irritation
Exposure time : 7 d
GLP : no

Acetic acid, potassium salt (1:1):

Species : Rabbit
Result : No eye irritation
Method : OECD Test Guideline 405

Respiratory or skin sensitization
Skin sensitization

Not classified due to lack of data.

Respiratory sensitization

Not classified due to lack of data.

Components:
Methane, 1,1'-sulfinylbis-:

Test Type : Local lymph node assay (LLNA)
Species : Mouse
Assessment : Does not cause skin sensitization.
Method : OECD Test Guideline 429
GLP : No information available.

Assessment : Mild eye irritation, Mild skin irritation

1,2,3-Propanetriol:

Assessment : Mild eye irritant, Mild respiratory irritant, No skin irritation

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020
Germ cell mutagenicity

Not classified due to lack of data.

Components:
Methane, 1,1'-sulfinylbis-

Genotoxicity in vitro : Test Type: Microbial mutagenesis assay (Ames test)
 Test system: Salmonella typhimurium
 Metabolic activation: with and without metabolic activation
 Method: OECD Test Guideline 471
 Result: negative
 GLP: No information available.

Test Type: Chromosome aberration test in vitro
 Test system: Chinese hamster ovary cells
 Metabolic activation: with and without metabolic activation
 Method: OECD Test Guideline 473
 Result: negative
 GLP: No information available.

Test Type: sister chromatid exchange assay
 Test system: Chinese hamster ovary cells
 Metabolic activation: with and without metabolic activation
 Method: OECD Test Guideline 479
 Result: negative
 GLP: No information available.

Genotoxicity in vivo : Test Type: In vivo micronucleus test
 Species: Rat (male and female)
 Cell type: Bone marrow
 Application Route: Intraperitoneal injection
 Dose: 200, 1000, 5000 mg/kg/d
 Method: OECD Test Guideline 474
 GLP: yes

1,2,3-Propanetriol:

Genotoxicity in vitro : Test Type: Microbial mutagenesis assay (Ames test)
 Test system: Salmonella typhimurium
 Metabolic activation: with and without metabolic activation
 Result: negative
 GLP: No information available.

Test Type: In vitro mammalian cell gene mutation test
 Test system: Chinese hamster ovary cells
 Metabolic activation: with and without metabolic activation
 Method: OECD Test Guideline 476
 Result: negative
 GLP: No information available.

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020**Carcinogenicity**

Not classified due to lack of data.

Components:**1,2,3-Propanetriol:**

Species : Rat, male and female
 Application Route : Oral
 Exposure time : 2 Years
 GLP : No information available.
 Remarks : No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

Reproductive toxicity

Not classified due to lack of data.

Components:**Methane, 1,1'-sulfinylbis-:**

Effects on fertility : Species: Rat, male and female
 Application Route: Oral
 Dose: 100, 300, 1000 mg/kg bw/day
 Fertility: NOAEL: 1,000 mg/kg body weight
 Method: OECD Test Guideline 421
 GLP: yes

Effects on fetal development : Species: Rat, female
 Application Route: Oral
 Dose: 200, 1000, 5000 milligram per kilogram
 Duration of Single Treatment: 10 d
 Developmental Toxicity: NOAEL: 1,000 mg/kg body weight
 Method: OECD Test Guideline 414
 GLP: yes

1,2,3-Propanetriol:

Effects on fertility : Test Type: Two-generation study
 Species: Rat, male and female
 Application Route: Oral
 Dose: 2000 mg/kg bw/day
 Fertility: NOAEL: 2,000 mg/kg body weight

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

GLP: no

Effects on fetal development : Species: Rabbit, female
 Application Route: Oral
 Dose: 11.8, 54.8, 254.5, 1180 mg/kg bw/day
 Duration of Single Treatment: 29 d
 Developmental Toxicity: NOAEL: 1,180 mg/kg bw/day
 GLP: no

STOT-single exposure

Not classified due to lack of data.

Components:**DNA-dependent DNA polymerase:**

Assessment : The substance or mixture is not classified as specific target organ toxicant, single exposure.

STOT-repeated exposure

Not classified due to lack of data.

Components:**DNA-dependent DNA polymerase:**

Assessment : The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Repeated dose toxicity**Components:****Methane, 1,1'-sulfinylbis-:**

Species : Monkey, male and female
 NOAEL : 2970 mg/kg
 NOAEL : 2,970 mg/kg
 Application Route : Oral
 Exposure time : 87 Weeks
 Dose : 990, 2970, 8910 mg/kg
 Method : OECD Test Guideline 452
 GLP : no

Species : Rat, male and female
 NOAEC : 2783 mg/l
 Application Route : Inhalation
 Test atmosphere : vapor
 Exposure time : 13 Weeks
 Dose : 0.310, 0.964, 2.783 mg/l
 Method : OECD Test Guideline 413
 GLP : yes

Species : Monkey, male and female
 NOAEL : > 8910 mg/kg
 NOAEL : > 8,910 mg/kg
 Application Route : Dermal

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 12. ECOLOGICAL INFORMATION
Ecotoxicity**Components:****Methane, 1,1'-sulfinylbis-:**

- Toxicity to fish : LC50 (Danio rerio (zebra fish)): > 25,000 mg/l
End point: mortality
Exposure time: 96 h
Test Type: static test
Analytical monitoring: yes
Method: OECD Test Guideline 203
GLP: yes
- Toxicity to daphnia and other aquatic invertebrates : EC50 (Daphnia magna (Water flea)): 24,600 mg/l
Exposure time: 48 h
Test Type: static test
Analytical monitoring: yes
Method: OECD Test Guideline 202
GLP: No information available.
- Toxicity to algae/aquatic plants : EC50 (Pseudokirchneriella subcapitata (green algae)): 17,000 mg/l
End point: Growth rate
Exposure time: 72 h
Test Type: static test
Analytical monitoring: yes
Method: OECD Test Guideline 201
GLP: yes
- Toxicity to microorganisms : EC50 (activated sludge): 10 - 100 mg/l
Exposure time: 0.5 h
Analytical monitoring: no
Method: ISO 8192
GLP: No information available.

Ecotoxicology Assessment

- Toxicity Data on Soil : Not expected to adsorb on soil.
- Other organisms relevant to the environment : No data available

1,2,3-Propanetriol:

- Toxicity to fish : LC50 (Oncorhynchus mykiss (rainbow trout)): 54,000 mg/l
End point: mortality
Exposure time: 96 h
Test Type: static test
GLP: no
- Toxicity to daphnia and other aquatic invertebrates : LC50 (Daphnia magna (Water flea)): 1,955 mg/l
End point: mortality

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020Exposure time: 48 h
Test Type: static test
Analytical monitoring: no
GLP: no

Toxicity to algae/aquatic plants : (Scenedesmus quadricauda (Green algae)): > 10,000 mg/l
End point: Growth rate
Exposure time: 8 d
Test Type: static test
GLP: no

Toxicity to microorganisms : EC50 (Pseudomonas putida): > 10,000 mg/l
End point: Growth rate
Exposure time: 16 h
Test Type: static test
GLP: No information available.

Ecotoxicology Assessment

Acute aquatic toxicity : This product has no known ecotoxicological effects.

Chronic aquatic toxicity : This product has no known ecotoxicological effects.

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Acetic acid, potassium salt (1:1):

Toxicity to fish : LC50 (Brachydanio rerio (zebrafish)): 496.35 mg/l
Exposure time: 96 h
Method: OECD Test Guideline 203

Toxicity to daphnia and other aquatic invertebrates : EC50 (Daphnia magna (Water flea)): 459.5 mg/l
Exposure time: 48 h
Method: OECD Test Guideline 202

Ecotoxicology Assessment

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

DNA-dependent DNA polymerase:
Ecotoxicology Assessment

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020
Persistence and degradability
Components:
Methane, 1,1'-sulfinylbis-:

Biodegradability : aerobic
 Inoculum: activated sludge
 Concentration: 2 mg/l
 Result: Not readily biodegradable.
 Biodegradation: 31 %
 Exposure time: 28 d
 Method: OECD Test Guideline 301D
 GLP: yes

1,2,3-Propanetriol:

Biodegradability : aerobic
 Inoculum: activated sludge
 Concentration: 226 mg/l
 Result: Readily biodegradable.
 Biodegradation: 94 %
 Exposure time: 24 h
 GLP: no

Acetic acid, potassium salt (1:1):

Biodegradability : Result: Readily biodegradable.

Bioaccumulative potential
Components:
Methane, 1,1'-sulfinylbis-:

Partition coefficient: n-
 octanol/water : log Pow: -1.35 (68 °F / 20 °C)
 pH: 7
 GLP: No information available.

1,2,3-Propanetriol:

Partition coefficient: n-
 octanol/water : log Pow: -1.75 (77 °F / 25 °C)
 pH: 7.4
 Method: OECD Test Guideline 107
 GLP: no

Acetic acid, potassium salt (1:1):

Bioaccumulation : Bioconcentration factor (BCF): 3.16
 Remarks: Due to the distribution coefficient n-octanol/water,
 accumulation in organisms is not expected.

Partition coefficient: n-
 octanol/water : Remarks: No data available

DNA-dependent DNA polymerase:

Partition coefficient: n- : Remarks: No data available



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

octanol/water

Mobility in soil

No data available

Other adverse effects

Product:

Ozone-Depletion Potential : Regulation: 40 CFR Protection of Environment; Part 82
Protection of Stratospheric Ozone - CAA Section 602 Class I
Substances
Remarks: This product neither contains, nor was
manufactured with a Class I or Class II ODS as defined by the
U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A +
B).

Additional ecological information : No data available

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : Can be disposed as waste water, when in compliance with
local regulations.

Contaminated packaging : Empty containers should be taken to an approved waste
handling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

Domestic regulation

49 CFR

Not regulated as a dangerous good

Special precautions for user

Remarks : Not dangerous goods in the meaning of ADR/RID, ADN,
IMDG-Code, ICAO/IATA-DGR

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 15. REGULATORY INFORMATION
CERCLA Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 304 Extremely Hazardous Substances Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

This material does not contain any components with a section 302 EHS TPQ.

SARA 311/312 Hazards : No SARA Hazards**SARA 313** : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.**Clean Air Act**

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

The following chemical(s) are listed under the U.S. Clean Air Act Section 111 SOCM I Intermediate or Final VOC's (40 CFR 60.489):

Methane, 1,1'-sulfinylbis-	67-68-5	>= 10 - < 20 %
1,2,3-Propanetriol	56-81-5	>= 10 - < 20 %

Clean Water Act

The following Hazardous Substances are listed under the U.S. CleanWater Act, Section 311, Table 116.4A:

Potassium hydroxide (K(OH))	1310-58-3	>= 0 - < 0.1 %
Glycine, N,N'-1,2-ethanediybis[N-(carboxymethyl)-	60-00-4	>= 0 - < 0.1 %

The following Hazardous Chemicals are listed under the U.S. CleanWater Act, Section 311, Table 117.3:

Potassium hydroxide (K(OH))	1310-58-3	>= 0 - < 0.1 %
Glycine, N,N'-1,2-ethanediybis[N-(carboxymethyl)-	60-00-4	>= 0 - < 0.1 %

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

This product does not contain any priority pollutants related to the U.S. Clean Water Act

US State Regulations**Massachusetts Right To Know**

1,2,3-Propanetriol	56-81-5
Sodium azide (Na(N3))	26628-22-8



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Pennsylvania Right To Know

Water	7732-18-5
Methane, 1,1'-sulfinylbis-	67-68-5
1,2,3-Propanetriol	56-81-5
Acetic acid, potassium salt (1:1)	127-08-2
Glycine, N-[2-hydroxy-1,1-bis(hydroxymethyl)ethyl]-	5704-04-1
Sodium azide (Na(N3))	26628-22-8
Potassium hydroxide (K(OH))	1310-58-3

Maine Chemicals of High Concern

Product does not contain any listed chemicals

Vermont Chemicals of High Concern

Product does not contain any listed chemicals

Washington Chemicals of High Concern

Product does not contain any listed chemicals

California Permissible Exposure Limits for Chemical Contaminants

1,2,3-Propanetriol	56-81-5
--------------------	---------

The ingredients of this product are reported in the following inventories:

- AIIC : Not in compliance with the inventory
- DSL : This product contains the following components that are not on the Canadian DSL nor NDSL.
 - NOTUPDATED_Uridine 5'-(tetrahydrogen triphosphate), 2'-deoxy-, disodium salt
 - NOTUPDATED_2'-Deoxyguanosine 5'-triphosphate sodium salt
 - 2'-Deoxycytidine 5'-triphosphate disodium salt
 - Primer / Oligonucleotide / Probe
 - Uridine, 2'-deoxy-, 5'-(trihydrogen pyrophosphate)
 - 2'-Deoxyguanosine 5'-(trihydrogen diphosphate)
 - NOTUPDATED_2'-Deoxyadenosine-5'-diphosphate
 - DNA-dependent DNA polymerase
- NZIoC : On the inventory, or in compliance with the inventory
- ENCS : Not in compliance with the inventory
- ISHL : Not in compliance with the inventory
- KECI : Not in compliance with the inventory
- PICCS : Not in compliance with the inventory

cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
 Date of first issue: 02-26-2020

- IECSC : Not in compliance with the inventory
- TCSI : Not in compliance with the inventory
- TSCA : Product contains substance(s) not listed on TSCA inventory.
- TECI : Not in compliance with the inventory

TSCA list


No substances are subject to a Significant New Use Rule.

No substances are subject to TSCA 12(b) export notification requirements.

SECTION 16. OTHER INFORMATION

Further information

NFPA 704:



HMIS® IV:

HEALTH	/	0
FLAMMABILITY	0	
PHYSICAL HAZARD	0	

HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. The "*" represents a chronic hazard, while the "/" represents the absence of a chronic hazard.

Full text of other abbreviations

- OSHA P0 : USA. Table Z-1-A Limits for Air Contaminants (1989 vacated values)
- OSHA Z-1 : USA. Occupational Exposure Limits (OSHA) - Table Z-1 Limits for Air Contaminants
- US WEEL : USA. Workplace Environmental Exposure Levels (WEEL)
- OSHA P0 / TWA : 8-hour time weighted average
- OSHA Z-1 / TWA : 8-hour time weighted average
- US WEEL / TWA : 8-hr TWA

AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response,

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 08-20-2025

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

US / Z8 / 2404



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 1. IDENTIFICATION

Product name : MMX-R1

Manufacturer or supplier's details

Company name of supplier : Roche Diagnostics Deutschland GmbH

Address : 116 Sandhoferstrasse
Mannheim, 68305
Deutschland

Telephone : +496217590

Telefax : +496217592890

E-mail address : info.dia-sds@roche.com

Emergency telephone

Im Notfall: : Werkschutzzentrale Roche +49(0)621-759-2203
Diagnostics GmbH

Giftnotruf: : Mainz +49(0)6131-19240

Recommended use of the chemical and restrictions on use

Recommended use : Laboratory chemicals
Refer to product literature for further details.

SECTION 2. HAZARDS IDENTIFICATION

GHS classification in accordance with the OSHA Hazard Communication Standard (29 CFR 1910.1200)

Not a hazardous substance or mixture.

GHS label elements

No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Other hazards

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Components

No hazardous ingredients

SECTION 4. FIRST AID MEASURES

General advice : Do not leave the victim unattended.



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- If inhaled : If unconscious, place in recovery position and seek medical advice.
If symptoms persist, call a physician.
- In case of skin contact : Wash off with soap and water.
- In case of eye contact : Remove contact lenses.
Protect unharmed eye.
If eye irritation persists, consult a specialist.
- If swallowed : Keep respiratory tract clear.
Do not give milk or alcoholic beverages.
Never give anything by mouth to an unconscious person.
If symptoms persist, call a physician.
- Most important symptoms and effects, both acute and delayed : None known.
- Notes to physician : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

- Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Unsuitable extinguishing media : High volume water jet
- Specific hazards during fire fighting : No information available.
- Hazardous combustion products : No hazardous combustion products are known
- Further information : Standard procedure for chemical fires.
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
- Special protective equipment for fire-fighters : Wear self-contained breathing apparatus for firefighting if necessary.

SECTION 6. ACCIDENTAL RELEASE MEASURES

- Personal precautions, protective equipment and : Refer to protective measures listed in sections 7 and 8.



cobas SARS-CoV-2 Qualitative

Version 4.0 Revision Date: 08-20-2025 Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

emergency procedures

Environmental precautions : No special environmental precautions required.

Methods and materials for containment and cleaning up : Wipe up with absorbent material (e.g. cloth, fleece).
Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

Advice on protection against fire and explosion : Normal measures for preventive fire protection.

Advice on safe handling : For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.

Conditions for safe storage : Electrical installations / working materials must comply with the technological safety standards.

Further information on storage conditions : See label, package insert or internal guidelines

Materials to avoid : No materials to be especially mentioned.

Further information on storage stability : No decomposition if stored and applied as directed.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Contains no substances with occupational exposure limit values.

Engineering measures : No data available

Personal protective equipment

Respiratory protection : No personal respiratory protective equipment normally required.

Hand protection

In case of contact through splashing:
Material : Nitrile rubber
Break through time : > 30 min
Glove thickness : > 0.11 mm

In case of full contact:
Material : butyl-rubber
Break through time : > 480 min
Glove thickness : > 0.4 mm



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

- Remarks : Wear appropriate protective gloves to prevent skin contact.
Replace torn or punctured gloves promptly.
- Eye protection : Safety glasses
- Skin and body protection : Protective suit
- Hygiene measures : General industrial hygiene practice.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

- Appearance : liquid
- Color : No data available
- Odor : No data available
- Odor Threshold : No data available
- pH : 6.4
- Melting point/ range : No data available
- Boiling point/boiling range : No data available
- Flash point : does not flash
- Evaporation rate : No data available
- Flammability (solid, gas) : The product is not flammable., Does not sustain combustion.
- Flammability (liquids) : Does not sustain combustion.
- Self-ignition : No data available
- Upper explosion limit / Upper flammability limit : No data available
- Lower explosion limit / Lower flammability limit : No data available
- Vapor pressure : No data available
- Relative vapor density : No data available
- Relative density : No data available
- Solubility(ies)
Water solubility : completely miscible



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Solubility in other solvents	:	No data available
Partition coefficient: n-octanol/water	:	No data available
Autoignition temperature	:	No data available
Decomposition temperature	:	No data available
Viscosity		
Viscosity, dynamic	:	No data available
Viscosity, kinematic	:	No data available
Explosive properties	:	Not explosive
Oxidizing properties	:	The substance or mixture is not classified as oxidizing.
Particle characteristics		
Particle Size Distribution	:	Not applicable

SECTION 10. STABILITY AND REACTIVITY

Reactivity	:	No dangerous reaction known under conditions of normal use.
Chemical stability	:	Stable under normal conditions.
Possibility of hazardous reactions	:	Stable under recommended storage conditions. No hazards to be specially mentioned.
Conditions to avoid	:	No data available
Incompatible materials	:	Not applicable
Hazardous decomposition products	:	No hazardous decomposition products are known.

SECTION 11. TOXICOLOGICAL INFORMATION

Acute toxicity

Not classified due to lack of data.

Skin corrosion/irritation

Not classified due to lack of data.

Serious eye damage/eye irritation

Not classified due to lack of data.



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Respiratory or skin sensitization

Skin sensitization

Not classified due to lack of data.

Respiratory sensitization

Not classified due to lack of data.

Germ cell mutagenicity

Not classified due to lack of data.

Carcinogenicity

Not classified due to lack of data.

IARC No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

Reproductive toxicity

Not classified due to lack of data.

STOT-single exposure

Not classified due to lack of data.

STOT-repeated exposure

Not classified due to lack of data.

Aspiration toxicity

Not classified due to lack of data.

Further information

Product:

Remarks : No data available

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Product:

Ecotoxicology Assessment

Toxicity Data on Soil : Not expected to adsorb on soil.

Other organisms relevant to the environment : No data available

Persistence and degradability

No data available



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

Product:

Ozone-Depletion Potential : Regulation: 40 CFR Protection of Environment; Part 82 Protection of Stratospheric Ozone - CAA Section 602 Class I Substances
Remarks: This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

Additional ecological information : No data available

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : Can be disposed as waste water, when in compliance with local regulations.
Contaminated packaging : Empty containers should be taken to an approved waste handling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

Domestic regulation

49 CFR

Not regulated as a dangerous good

Special precautions for user

Remarks : Not dangerous goods in the meaning of ADR/RID, ADN, IMDG-Code, ICAO/IATA-DGR

cobas SARS-CoV-2 Qualitative
Version
4.0Revision Date:
08-20-2025Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

SECTION 15. REGULATORY INFORMATION
CERCLA Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 304 Extremely Hazardous Substances Reportable Quantity

Listed substances in the product are at low enough levels to not be expected to exceed the RQ

SARA 302 Extremely Hazardous Substances Threshold Planning Quantity

This material does not contain any components with a section 302 EHS TPQ.

SARA 311/312 Hazards : No SARA Hazards**SARA 313** : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.**Clean Air Act**

This product neither contains, nor was manufactured with a Class I or Class II ODS as defined by the U.S. Clean Air Act Section 602 (40 CFR 82, Subpt. A, App.A + B).

This product does not contain any hazardous air pollutants (HAP), as defined by the U.S. Clean Air Act Section 112 (40 CFR 61).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 112(r) for Accidental Release Prevention (40 CFR 68.130, Subpart F).

This product does not contain any chemicals listed under the U.S. Clean Air Act Section 111 SOCM Intermediate or Final VOC's (40 CFR 60.489).

Clean Water Act

The following Hazardous Substances are listed under the U.S. CleanWater Act, Section 311, Table 116.4A:

Acetic acid	64-19-7	>= 0 - < 0.1 %
Potassium hydroxide (K(OH))	1310-58-3	>= 0 - < 0.1 %

The following Hazardous Chemicals are listed under the U.S. CleanWater Act, Section 311, Table 117.3:

Acetic acid	64-19-7	>= 0 - < 0.1 %
Potassium hydroxide (K(OH))	1310-58-3	>= 0 - < 0.1 %

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

This product does not contain any priority pollutants related to the U.S. Clean Water Act

US State Regulations**Massachusetts Right To Know**

Sodium azide (Na(N3)) 26628-22-8

Pennsylvania Right To Know

Water	7732-18-5
Acetic acid, manganese(2+) salt, tetrahydrate	6156-78-1
Sodium azide (Na(N3))	26628-22-8

Maine Chemicals of High Concern

Product does not contain any listed chemicals



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Vermont Chemicals of High Concern

Product does not contain any listed chemicals

Washington Chemicals of High Concern

Product does not contain any listed chemicals

The ingredients of this product are reported in the following inventories:

- AIIC : On the inventory, or in compliance with the inventory
- DSL : This product contains the following components that are not on the Canadian DSL nor NDSL.

Acetic acid, manganese(2+) salt, tetrahydrate
- NZIoC : On the inventory, or in compliance with the inventory
- ENCS : Not in compliance with the inventory
- ISHL : Not in compliance with the inventory
- KECI : On the inventory, or in compliance with the inventory
- PICCS : On the inventory, or in compliance with the inventory
- IECSC : On the inventory, or in compliance with the inventory
- TCSI : On the inventory, or in compliance with the inventory
- TSCA : Product contains substance(s) not listed on TSCA inventory.
- TECI : Not in compliance with the inventory

TSCA list

No substances are subject to a Significant New Use Rule.

No substances are subject to TSCA 12(b) export notification requirements.

SECTION 16. OTHER INFORMATION

Further information


cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

NFPA 704:



HMIS® IV:

HEALTH	/	0
FLAMMABILITY		0
PHYSICAL HAZARD		0

HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. The "*" represents a chronic hazard, while the "/" represents the absence of a chronic hazard.

Full text of other abbreviations

AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CERCLA - Comprehensive Environmental Response, Compensation, and Liability Act; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DOT - Department of Transportation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; EHS - Extremely Hazardous Substance; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; HMIS - Hazardous Materials Identification System; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; MSHA - Mine Safety and Health Administration; n.o.s. - Not Otherwise Specified; NFPA - National Fire Protection Association; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; RCRA - Resource Conservation and Recovery Act; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RQ - Reportable Quantity; SADT - Self-Accelerating Decomposition Temperature; SARA - Superfund Amendments and Reauthorization Act; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECl - Thailand



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Revision Date : 08-20-2025

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

US / Z8 / 2404

SAFETY DATA SHEET



cobas SARS-CoV-2 Qualitative

Version
4.0

Revision Date:
08-20-2025

Date of last issue: 08-20-2025
Date of first issue: 02-26-2020

