



a member of the Roche group

## **CERTIFICATE OF ANALYSIS**

---

Product Name: NeuVision™ Blocking Reagent  
Ventana REF: 760-1028  
Roche Number: 05266220001  
Batch Number: D02070

Test	Acceptance Criteria	Result
Functional Staining on breast carcinoma tissue and reference tissue slide	Reference and Test Slide- acceptable staining with minimal background.	PASS
	Non-Blocking Control Slide- equivalent or less background than slides stained without blocking reagent.	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of intrinsic particulate matter.	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL.	PASS

---

This lot met the specifications.

Expiry Date: 2014-04-28

Manufacture Date: 18-APR-2013

Angela Glenn  22-Apr-2013  
Quality Department Date