



a member of the Roche group

CERTIFICATE OF ANALYSIS

Product Name: VENTANA Basal Cell Cocktail (34βE12+p63)

Ventana REF: 790-4536 or 790-1010

Roche Number 06364497001 or 06419445001

Batch Number: G05632

Test	Acceptance Criteria	Result
Functional Staining on normal skin and prostate carcinoma	The test batch must score 3.0+ or greater on normal skin and normal prostate CA tissues	PASS
	The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result	PASS
	The background intensity must be 0.5+ or less	PASS
	The intensity of the negative reagent control slide must be 0.5 or less	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiration date: 2018-07-28

Manufacture date: 28-Jul-2016

Ken Eblen

AVG/2016

Quality Department

Date