



Certificate of Analysis

Product name

KIT COBAS 4800 CMV 120T CE-IVD

Material No.

07865970190

Lot No.

E18065

Date of Manufacturing

2018-09-24

Expiry Date

2019-09

Parameter	Result	Unit
Components	1	
5XLOD-Sens Reps-Pos	100	%
C48CMVKIT-120T	PASS	

The Visual Inspection should be as follows:

10 X ADAPTER c4800 CMV MMX R2, 0.5 mL
Black microtube and cap with a gold insert, fitted firmly into a labeled adapter.

10 X ADAPTER c4800 MMX R1, 1.75 mL
Natural microtube and cap with a teal insert, fitted firmly into a labeled adapter.

10 X ADAPTER c4800 DNA QS, 1.75 mL
Natural microtube and cap with a purple insert, fitted firmly into a labeled adapter.

Calibration Coefficient Verification
Kit Calibration Coefficient values encoded in the Assay

Issue date 2018-09-27

Page 1 of 2

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771



Certificate of Analysis

Product name

KIT COBAS 4800 CMV 120T CE-IVD

Material No.

07865970190

Lot No.

E18065

Date of Manufacturing

2018-09-24

Expiry Date

2019-09

Parameter**Result****Unit**

Value Card must match the values recorded in SAP.

CALHI1-Accuracy	100	%
CALHI1-Prec-SD	0.02	
STDLO-Accuracy	100	%
STDLO-Prec-SD	0.03	
NHP-Negativity-NR	100	%

This lot meets our specification

Disposition by

Mujib Walli

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2018-09-27

Page 2 of 2

Roche Molecular Systems, Inc.
Quality Management
1080 U.S. Highway 202 South
Branchburg, NJ 08876-3771