



Certificate of Analysis

Product Name: CD21 (EP3093) Rabbit Monoclonal Antibody
Catalog Number(s): 760-4438
Roche Part Number: 05973902001
Lot Number: V0000585
Date of Manufacture: 2016-10-12
Expiration Date: 2019-09-30

| Test | Acceptance Criteria | Results |
|---|--|--|
| Functional testing on the following tissue(s): Lymph Node, Tonsil | Intensity difference between sample and a reference must be less than or equal to 0.5 | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| | Background difference between sample and a reference must be less than or equal to 0.5 | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| Appearance | Transparent liquid and free of visible foreign material. | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| Bioburden | Visible bacterial colonies must be less than 100CFU/mL. | <input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL |
| ELISA | 0.01-1.0 µg/mL | 0.36 µg/mL |
| Source | N/A | Supernatant |

Additional information: N/A



CELL MARQUE
6600 Sierra College Blvd.
Rocklin, California 95677
(916) 746-8900 Fax: (916) 746-8989

Authorized Signature:  Date: 7/27/17
Quality Assurance Representative

Issued: 8/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]

Manufactured by: Cell Marque Corporation
6600 Sierra College Boulevard • Rocklin, California 95677 • USA • Telephone (916) 746-8900 • Fax (916) 746-8976