



Certificate of Quality

Product Name: 6-Acetylmorphine Enzyme Immunoassay
Product Size: c pack kit 240 tests (24.0 mL R₁*, 10.1 mL R₂*)
* overfill volume included

Reference Number: 09067795190

Description: The 6-Acetylmorphine Enzyme Immunoassay is a homogeneous enzyme immunoassay intended for use in the qualitative and semi-quantitative determination of 6-acetylmorphine in human urine.

Lot Number: 25034309



Date of Manufacture: October 27, 2025

Expiration Date: April 30, 2027

QC/QA results: cobas c 501 standard curve results: mA/min

Negative (0 ng/mL)	:	348
Low (5 ng/mL)	:	381
Negative Control (7.5 ng/mL)	:	394
Cutoff (10 ng/mL)	:	417
Positive Control (12.5 ng/mL)	:	432
Intermediate (20 ng/mL)	:	481
High (40 ng/mL)	:	588

Storage/Preparation: The reagents are ready to use. No reagent preparation is required. All assay components should be stored at 2-8°C when not in use.

 
QC/QA: Jason Nguyen and Navpreet Pal

Date: January 08, 2026