



Certificate of Analysis

Product name

VENTANA HE 600 Eosin

Material No.

06544304001

Lot No.

H16508

Date of Manufacturing

2021-06-24

Expiry Date

2023-06-23

Parameter**Result****Unit**

Final Acceptance Testing

Pass

Small intestine

Differential staining of red blood cells, smooth muscle, and collagen must be present.

The intensity of the Eosin staining in the test slide must be more intense than the low-intensity Eosin rated reference slides and less intense than the high-intensity Eosin rated reference slides.

It is unacceptable for the test slide to be equivalent in intensity to the low- or high-intensity rated reference slides.

At least 3 of the 4 small intestine sections on the test slide must meet criteria.

This lot meets our specification

Disposition by

Antonio Reyes

Issue date 2021-07-10

Page 1 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 520 887 2155
<http://www.ventana.com>



Certificate of Analysis

Product name

VENTANA HE 600 Eosin

Material No.

06544304001

Lot No.

H16508

Date of Manufacturing

2021-06-24

Expiry Date

2023-06-23

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2021-07-10

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 520 887 2155
<http://www.ventana.com>