

Certificate of Analysis

**Product name**

CONFIRM[®] anti-Topoisomerase II α (JS5B4) Rabbit Monoclonal Primary Antibody

Material No.

05479339001

Lot No.

G07622

Date of Manufacturing

28-Nov-2016

Expiry Date

2018-11-28

Parameter**Result****Unit**

Visual Inspection

Pass

Reagent must be free of foreign material and substantially free of particulate matter

Reagent must be clear and free of turbidity

Functional Staining

Pass

Breast Cancer

The test batch must score 3.5 or greater on breast cancer tissue

The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result

The background intensity must be within 0.5 of the negative control slide

The intensity of the negative reagent control slide must be 0.5 or less

Bioburden

Pass

Visible bacterial colonies must be 100 or less CFU/mL

This lot meets our specification

Disposition by

Lucia Romero

Lucia Romero LB

Issue date

Dec/2016