



# Certificate of Analysis

**Product name**

PROTEASE 1

**Material No.**

05266688001

**Lot No.**

87681801

**Date of Manufacturing**

2025-03-25

**Expiry Date**

2027-01-15

**Parameter****Result****Unit**

Visual Inspection

Corresponds

Product should be clear and colorless, but it is acceptable for the reagent to contain a white flocculent material of a non-biological origin.

Absorbance

Corresponds

Reference Sample and Test Sample Mean Absorbance Unit must be in the range of 0.50 and 1.30.

Bioburden

Corresponds

Visible bacterial colonies must be 100 or less CFU/mL

This lot meets our specification

**Disposition by**

Jie Wang

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2026-01-14

Page 1 of 1

Roche Diagnostics (Suzhou) Ltd.

No.259 Zhongyuan Road, 215028 Suzhou, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-0512-6999 8200; Fax: +86-0512-6999 8222;

罗氏诊断产品(苏州)有限公司

中国江苏苏州工业园区钟园路259号