



a member of the Roche group

## CERTIFICATE OF ANALYSIS

Product: Confirm™ anti- CD68 (KP-1) Primary Antibody  
Ventana REF: 790-2931  
Roche Number: 05278252001  
Batch Number: C07583

Test	Acceptance Criteria	Result
Functional Staining on Tonsil, Spleen, or Lymphoma	Reference slide reading must score 3.5+ or greater on a scale of 0-4+	PASS
	Background reading must score 0.5+ or less on a scale of 0-4+	PASS
	Intensity difference between sample and a reference must be less than or equal to 0.5+	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of intrinsic particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiry date: 2014-03-05

Manufacture date : 21-Sep-2012

Angela Glenn  31 OCT 2012  
Quality Department Date