



## Certificate of Analysis

**Product Name:** IgM Rabbit Polyclonal Antibody  
**Catalog Number(s):** 760-2654  
**Roche Part Number:** 05267641001  
**Lot Number:** V0000103  
**Date of Manufacture:** 2016-07-28  
**Expiration Date:** 2019-06-30

Test	Acceptance Criteria	Results
Functional testing on the following tissue(s): Tonsil	Intensity difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
	Background difference between sample and a reference must be less than or equal to 0.5	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Appearance	Transparent liquid and free of visible foreign material.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
Bioburden	Visible bacterial colonies must be less than 100CFU/mL.	<input checked="" type="checkbox"/> PASS <input type="checkbox"/> FAIL
ELISA	1.0-10.0 µg/mL	3.28 µg/mL
Source	N/A	Anti-sera

Additional information: N/A



Authorized Signature:  Date: 10/17/17  
Quality Assurance Representative

Issued: 6/3/16 Rev 1.0

Document #QF0901-008-02 [CM Legacy #]

Manufactured by: Cell Marque Corporation  
6600 Sierra College Boulevard • Rocklin, California 95677 • USA • Telephone (916) 746-8900 • Fax (916) 746-8976