

**CERTIFICATE OF ANALYSIS**


---

Product: Confirm™ anti-Progesterone Receptor (PR)(1E2) Rabbit Monoclonal Primary Antibody  
 Ventana REF: 790-2223  
 Roche Number: 05277990001  
 Batch Number: E07088

| Parameter   | Result  |      |
|---|---|------|
| Functional Staining on breast carcinoma tissue and reference tissue slide | Nuclear staining must be present in 51% or more of the tumor cells                              | PASS |
|   | Background staining reading must score 0.5+ or less on a scale of 0-4+                          | PASS |
|   | Intensity difference between sample and a reference must be less than or equal to 0.5           | PASS |
| Visual Inspection   | Reagent must be free of foreign material and substantially free of intrinsic particulate matter | PASS |
|   | Reagent must be clear and free of turbidity   | PASS |
| Bioburden   | Visible bacterial colonies must be less than or equal to 100 CFU/mL                             | PASS |

This lot met the specifications.

Expiry date: 2016-08-28

Manufacture date: 23-Sep-2014

Maria Delia Lamadrid *Maria Delia Lamadrid*

14-OCT-2014

**Quality Department**

**Date**