



Certificate of Analysis

Product name
PATHWAY anti-HER2/neu (4B5) Rabbit Monoclonal Primary Antibody

Material No.
05278368001

Lot No.	Date of Manufacturing	Expiry Date
H00315	15-JAN-2021	2022/07/11

Parameter	Result	Unit
Visual Inspection	Pass	

Reagent must be free of foreign material and substantially free of particulate matter

Reagent must be clear and free of turbidity

Functional Staining	Pass	
---------------------	------	--

Breast Carcinoma Tissue, Gastric Carcinoma Tissue, and 4 in 1 Cell Line.

The specific staining of the test batch must be acceptable on all tissues and cell line cores using both the *View* DAB IHC and the *ultraView* DAB IHC detection kits.

The test batch background intensity must not exceed 0.5 pt on any tissue or cell line core.

All tissues and cell line cores stained with the negative control reagent must be negative for specific staining

Bioburden	Pass	
-----------	------	--

Visible bacterial colonies must be 100 or less CFU/mL

This lot meets our specification

Disposition by

 NICOLE CAREY

Issue date

15 FEB 2021