

CERTIFICATE OF ANALYSIS

Product Name: Kappa DNP Probe (637-666), ASR
 Ventana REF: 760-1204
 Roche Number: 06543421001
 Batch Number: D03254


Test	Acceptance Criteria	Result
Functional Staining on qualified normal tonsil tissue slide	Moderate blue cytoplasmic staining in plasma cells	PASS
	Blue staining of the Kappa DNA spot	PASS
	Negative staining of the Lambda DNA spot	PASS
	Background must not interfere with signal reading of plasma cells	PASS
Visual Inspection	Reagent must be clear and free of turbidity	PASS

This lot met the specifications.

Expiry Date: 2014-05-28

Manufacture Date: 09-MAY-2013

Elizabeth Garcia



17-May-2013

Quality Department

Date