



Certificate of Analysis

Product name

CONFIRM ANTI-KI-67 (30-9) RABBIT MONOCLO

Material No.

05278384001

Lot No.

F07946

Date of Manufacturing

2019-05-06

Expiry Date

2021-05-05

Parameter**Result****Unit**

Final Acceptance Testing

Pass

Breast carcinoma, normal tonsil, or normal lymph node.

The test batch must stain 51% or greater of cells on breast carcinoma, normal tonsil, or normal lymph node tissue.

The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result.

The intensity of the negative reagent control slide must be 0.5 or less.

Bioburden

Pass

Visible bacterial colonies must be 100 or less CFU/mL

Visual Inspection

Pass

Reagent must be clear, and free of turbidity.

Reagent must be free of foreign material and

Issue date 2019-06-04

Page 1 of 2



Certificate of Analysis

Product name

CONFIRM ANTI-KI-67 (30-9) RABBIT MONOCLO

Material No.

05278384001

Lot No.

F07946

Date of Manufacturing

2019-05-06

Expiry Date

2021-05-05

Parameter**Result****Unit**

substantially free of intrinsic particulate matter.

This lot meets our specification

Disposition by

Ken Eblen

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2019-06-04

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 520 887 2155
<http://www.ventana.com>