



a member of the Roche group

CERTIFICATE OF ANALYSIS

Product: GMS II Staining Kit
Ventana **REF**: 860-028
Roche Number: 05412749001
Batch Number: G06273

Test	Acceptance Criteria	Result
Functional Staining on Pneumocystis carinii and Aspergillus sp. fungus	Fungus should stain black against a light green background	PASS
	Aspergillus sp. Fungal hyphae should appear as black filaments	PASS
	Pneumocystis carinii should appear as clusters of black spots	PASS

This lot met the specifications.

Expiry Date: 2018-02-28

Manufacture Date : 28-AUG-2016

Ken Eblen 

SEP/2016

Quality Department

Date