



Certificate of Analysis

Product name

ISH AMPLIFICATION KIT

Material No.

07729022001

Lot No.

N03252

Date of Manufacturing

2025-03-24

Expiry Date

2026-08-20

Parameter**Result****Unit**

Visual Inspection

Pass

Reagent must be clear of foreign material and particulate matter, and free of turbidity

Functional Staining

Pass

ISH Amplification Kit test slides are read and compared to OptiView DAB Detection Kit test slides.

Specificity is acceptable for CONFIRM anti-bcl-2 tested with OptiView DAB Detection Kit. Interfollicular T cells in Tonsil, B cells of the Mantle Zone. Localization is cytoplasmic and when tested with ISH Amplification Kit is greater than or equal to 1.0 point when compared to OptiView DAB Detection Kit.

Background Grading is acceptable when no Antibody Control (Reaction Buffer) tested with OptiView DAB Detection Kit with some discernable non-specific staining, which would not be obtrusive to stain interpretation nor appear as a false positive signal. Tested with ISH Amplification Kit are acceptable when equivalent to slides tested with OptiView DAB Detection Kit.

This lot meets our specification

Issue date 2025-04-09

Page 1 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>



Certificate of Analysis

Product name

ISH AMPLIFICATION KIT

Material No.

07729022001

Lot No.

N03252

Date of Manufacturing

2025-03-24

Expiry Date

2026-08-20

Disposition by

Lupe Villarreal

This certificate is a computer printout and has therefore not been signed by hand.

Issue date 2025-04-09

Page 2 of 2

Ventana Medical Systems, Inc.
1910 Innovation Park Drive
Tucson, Arizona 85755
Phone: + 1 800 227 2155
<http://www.ventana.com>