



# Certificate of Analysis

**Product name**

Creatine Kinase (CK), 100mg

**Material No.**

10127566001

**Lot No.**

89507828

**Date of Manufacturing**

2025-12-03

**Expiry Date**

2026-12

Parameter	Result	Unit
Appearance white to slightly pink	corresponds	
or yellowish lyo.	corresponds	
stabilized	corresponds	
Act.(25 deg.C;	412	u /mgL
with creatine-P a. ADP;NAC as activator)	corresponds	
Contaminating activities	corresponds	
ATPase	< 0,001	%
HK(with PK/LDH)	< 0,0005	%
Myokinase	< 0,001	%

This lot meets our specification

This certificate is a computer printout and has therefore not been signed by hand.

**Issue date** 2025-12-03

Page 1 of 1

Roche Diagnostics GmbH

Sandhofer Straße 116; D-68305 Mannheim; Telefon +49-621-759-0; Telefax +49-621-759-2890

Sitz der Gesellschaft: Mannheim - Registergericht: AG Mannheim HRB 3962 - Geschäftsführung: Dr. Claudia Fleischer;

Dr. Virginia Bastian - Aufsichtsratsvorsitzender: Dr. Thomas Schinecker