

**CERTIFICATE OF ANALYSIS**


---

Product Name: Anti-c-MYC (Y69) Rabbit Monoclonal Primary Antibody

Ventana [REF]: 790-4628

Roche Number: 06504612001

Batch Number: G01437

Test	Acceptance Criteria	Result
Functional Staining on lymphoma	The test batch must score 2.0+ or greater on lymphoma tissue	PASS
	The specific stain intensity must also be within 0.5 pt. of the reference batch for a passing result	PASS
	The background intensity must be within 0.5 of the negative control slide	PASS
	The intensity of the negative reagent control slide must be 0.5 or less	PASS
Visual Inspection	Reagent must be free of foreign material and substantially free of particulate matter	PASS
	Reagent must be clear and free of turbidity	PASS
Bioburden	Visible bacterial colonies must be less than or equal to 100 CFU/mL	PASS

This lot met the specifications.

Expiration date: 2018/03/28

Manufacture date: 28/MAR/2016

  
Quality Department

Jeremy Thomas

AUG 1/2016

Date